

NEXUS ARCH

A MULTICENTER STUDY EVALUATING THE INITIAL EXPERIENCE WITH A NOVEL AORTIC ARCH STENT GRAFT SYSTEM

Ann Surg 2021 – accepted

Published Feb 2023

PATIENT TREATED WITH NEXUS™ FOR AORTIC ARCH ANEURYSM*


Baseline



1-Month



12-Months



**Courtesy of Dr Andrew Hill, Auckland New Zealand.*

STUDY DESIGN



Prospective cohort
28 patients in **5** centers



From first-in-man
(FIM n=18) study
(NCT02365454) or as
compassionate use
(n=10) with systematic
data collection
(NCT03420066)

	First-in-human (n=18)	Compassionate use (n=10)
University Hospital Zürich, Switzerland	14	2
Azienda Complesso Ospedaliero San Filippo Neri, Rome, Italy	3	2
University Hospital Auckland, New Zealand	-	1
Toronto General Hospital, Canada	-	5
University Hospital, Düsseldorf, Germany	1	-


BASELINE CHARACTERISTICS

Variable	First in Man Cohort (n = 18)	Compassionate Cohort (n = 10)	Entire Cohort (n = 28)	P-value
Age (years), mean ± SD	71.7 ± 5.9	73.2 ± 6.9	72.2 ± 6.2	0.55
NEXUS procedure time (min) [†]	16/18 (88.9%)	6/10 (60%)	22/28 (78.6%)	0.15
BMI kg/m ² , mean ± SD	28.6 ± 5.4	27.2 ± 7.2	28.1 ± 6.0	0.57
DM	2/18 (11.1%)	0/10	2/28 (7.1%)	0.52
HTN	18/18 (100%)	9/10 (90%)	27/28 (96.4%)	0.36
DLP	6/18 (33.3%)	5/9 (55.5%)	11/27 (40.7%)	0.41
Current smoker	3/18 (16.7%)	1/10 (10%)	4/28 (14.3%)	1.00
COPD	5/18 (27.8%)	4/10 (40%)	9/28 (32.1%)	0.68
CAD	7/18 (38.9%)	2/10 (20%)	9/28 (32.1%)	0.42
Arrhythmia	6/18 (33.3%)	1/10 (10%)	7/28 (25%)	0.36
Previous sternotomy	12/18 (27.8%)	3/10 (30%)	15/28 (53.6%)	0.11
CVA/TIA	1/18 (5.6%)	1/10 (10%)	2/28 (7.1%)	1.00
CHF (NYHA III/IV)	0	2/10 (20%)	2/28 (7.1%)	0.12
PVD	2/18 (11.1%)	1/10 (10%)	3/28 (10.7%)	1.00
CKD (Creatine > 2.0)	2/18 (11.1%)	0	2/28 (7.1%)	0.52
Anemia	7/18 (38.9%)	2/10 (20%)	9/18 (32.1%)	0.42
ASA risk score ≥ 3	16/17 (94.1%)	9/10 (90%)	25/27 (92.5%)	0.70

ASA indicates American Society of Anesthesiology; CAD, coronary artery disease; CHF, congestive heart failure; CKD, chronic kidney disease; COPD, chronic obstructive lung disease; CVA, cerebrovascular accident; DLP, dyslipidemia; DM, diabetes mellitus; HTN, hypertension; IQR, interquartile range; NYHA, New York Heart Association; PVD, peripheral vascular disease; TIA, Transient ischemic attack.

CAUTION: Investigational Device – Limited by United States law to investigational use. Endospans devices bear the CE marking of conformity.

TREATED AORTIC ARCH ACCORDING TO MOST PROXIMAL PATHOLOGY*



*Roselli et al. Zone zero thoracic endovascular aortic repair:
A proposed modification to the classification of landing zones.
J Thorac Cardiovasc Surg. 2018;155(4):1381-9

PROCEDURAL CHARACTERISTICS

Variable	Median [IQR]
Total procedure time (min)*	185 [148–254]
NEXUS procedure time (min)†	80 [46.5–113]
Fluoroscopy time (min)‡	48 [37.5–54]
Contrast volume (mL)§	122.5 [102.5–187.5]
Length of hospitalization following procedure (d)	8.5 [7.0–14.7]
Number of patients admitted to ICU	16/28 (57.1%)
Length of ICU stay (d)	1.0 [1.0–3.0]
Complete percutaneous access	19/28 (67.9%)

*Total procedure time: skin-to-skin; data was collected for 26 patients

†NEXUS procedure time: time from access with the arch stent graft delivery system until retrieval of the ascending stent graft delivery system; data was collected for 17 patients

‡Fluoroscopy time data was collected for 25 patients

§Contrast volume data was collected for 24 patients

ICU indicates intensive care unit

MORTALITY & STROKE RATE



	30 Days	1 Year
Mortality Rate	7.1%	10.7%¹
Stroke		
<i>Disabling</i>	0.0%	0.0%
<i>Non-Disabling</i>	3.6%	7.1%²

1

Not device or aneurysm-related death per CEC Adjudication

2

Reported Stroke or TIA that completely recovered



Kaplan-Meier survival curves

A: Freedom from device-related unplanned reinterventions

B: Overall survival


ENDOLEAKS WITHOUT CLINICAL SEQUELAE

Type Ic
None

Type III

1 patient at 30d (resolved spontaneously) and 1 patient at 6 months – stable at 1 year – no aneurysm growth

Type Ia
None



Type II
4 non occluded LSA (expected) – successful staged embolization

Type Ib
1 patient – expected – pre-planned distal thoracic stent graft extension

OTHER OUTCOMES

A 1-year follow-up CTA

No graft migration, stent graft separation, branch occlusion, stent fracture, graft infoldings or collapse.

1 report of asymptomatic surgical graft occlusion (LCC artery to LSA)

No reports of periscope occlusion

Periscope stent grafts in 10 patients

Unplanned periscope interventions:

1 at 30 days and **2 at 12 months**

(Aneurysm enlargements >5mm)

Planned device related interventions: (periscope):

1 at 12 months

No stent migration

CONCLUSION

1

In patients with aortic arch pathologies that require landing in the ascending aorta, endovascular repair with the NEXUS Aortic Arch Stent Graft System can be performed with high success rate and promising results at 1 year.

2

Further follow-up is required to establish the long-term safety and effectiveness of this device.

28 Patients with 1 year FU

KEY OUTCOMES

100% procedural success.

At one month, the vascular pathology was excluded in **25 of 26** alive patients (**96.1%**). One year mortality was **10.7%**, without device or aneurysm related death.

Two patients (7.1%) reported non-disabling stroke or TIA at one year that recovered completely.

Three patients (10.7%) had device related unplanned reinterventions through one year.

NEXUS™ VS TERUMO

	NEXUS™: OUS (n=28) ^{1,2}	NEXUS™: US IDE (n=10) ^{1,3}	TERUMO (n=12)
Early Mortality	2/28 (7.1%)	1/10 (10%)	2/12 (17%)
Stroke	1/28 (3.6%)	0/10 (0%)	6/12 (50%)
Paralysis	0/28 (0%)	0/10 (0%)	1/12 (8%)
Respiratory Failure	1/28 (3.6%)	0/10 (0%)	1/12 (8%)

Outcomes

- Death: Acute 2 (17%), day 11 & 34
- Stroke: Late 1 (8%), day 246
- Stroke: All 6 (50%)
- Stroke: Disabling 1 (8%)
- Paralysis 1 (8%)
- MI 1 (8%)
- Respiratory Failure 1 (8%)
- Aortic Dissection 1 (8%)

LIMITATION: Arch Disease is High Risk Population – esp Aneurysm

Cleveland Clinic

Presented by Dr. Roselli
from the Cleveland Clinic at VEITH in 2019

¹ Results includes Bypass and NEXUS Implant

² NCT02365454; NCT03420066

³ NCT04471909