(en) ENGLISH

## Instructions for Use

ARTIVION™ | BioGlue®

### BioGlue Surgical Adhesive

[www.artivion.com/eifu/bioglue](http://www.artivion.com/eifu/bioglue)

A printed version of the Instructions for Use will be provided within seven days upon request to BioGlue customer service through any of the contact means listed below.

Phone: 888.427.9654 • Fax: 770.590.3753

Email: [customerservice.us@artivion.com](mailto:customerservice.us@artivion.com)

ARTIVION™

1655 Roberts Boulevard, NW  
Kennesaw, Georgia 30144, USA

Phone: +1 (770) 419-3355

FAX +1 (770) 590-3753

[www.artivion.com](http://www.artivion.com)

JOTEC GmbH

Lotzenäcker 23

72379 Hechingen, Germany

CE 0124

L09418.001 (2023-06-05)

## EXPLANATION OF SYMBOLS

|   |                                  |   |  |
|---|----------------------------------|---|--|
|    | Manufacturer                     |    | Consult instructions for use / Consult electronic instructions for use                     |
|    | Date of manufacture              |    | Caution  |
|    | EU Authorized Representative     | <b>RxONLY</b>   | Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician |
|    | Importer                         |    | Non-pyrogenic  |
|    | Medical Device                   |    | Contains biological material of animal origin  |
|    | Do not use if package is damaged |    | MR safe  |
|    | Single sterile barrier system    |    | Catalogue number   |
|    | Double sterile barrier system    |    | Batch code   |
|    | Sterilized using ethylene oxide  |    | Unique device identification   |
|  | Sterilized using irradiation     |  | Use-by date  |
|  | Do not re-sterilize              |  | Country of manufacture   |
|  | Do not re-use                    |  | Quantity   |
|  | Temperature limit                |   |  |

## PRODUCT DESCRIPTION

BioGlue® Surgical Adhesive is a two-component surgical adhesive composed of solutions of purified bovine serum albumin and glutaraldehyde. Once dispensed, the adhesive solution (in predefined ratio) are mixed within the applicator tip where cross-linking begins. The glutaraldehyde molecules covalently bond (cross-link) the BSA molecules to each other and, upon application, to the tissue proteins at the repair site, creating a flexible mechanical seal independent of the body's clotting mechanism. BioGlue Surgical Adhesive (hereinafter BioGlue) begins to polymerize within 20 to 30 seconds and reaches its bonding strength within 2 minutes. BioGlue also adheres to synthetic graft materials via mechanical interlocks within the interstices of the graft matrix.

The following accessories are sold separately to aid in the delivery of BioGlue Surgical Adhesive:

| Product Code | Product Description            |
|--------------|--------------------------------|
| BGAT-SY      | Applicator tip                 |
| BGAT-10-SY   | Syringe extender tip - 10cm    |
| BGAT-27-SY   | Syringe extender tip - 27cm    |
| BGST-12      | Spreader applicator tip - 12mm |
| BGST-16      | Spreader applicator tip - 16mm |
| BGDTE-10     | Delivery tip extension - 10cm  |
| BGDTE-27     | Delivery tip extension - 27cm  |
| BGDTE-35     | Delivery tip extension - 35cm  |

BioGlue syringes are available in 3 configurations - 2mL, 5mL, and 10mL. Each syringe is composed of bovine serum albumin (BSA) and glutaraldehyde solutions in a 4:1 ratio, respectively. The BSA solution is amber in color and free-flowing. The glutaraldehyde solution is clear and also free-flowing.

The specification for the BSA solution is a 45% (weight/volume ratio) solution. The maximum 45% BSA solution target weights for each size are: 2.71 grams (2mL syringe), 4.75 grams (5mL syringe), and 9.50 grams (10mL syringe). Based on these targets, the maximum quantity of animal origin material coming into contact with the patient when using a single device is 1.22 grams (2mL syringe), 2.14 grams (5mL syringe), and 4.23 grams (10mL syringe) for each configuration.

The specification for the glutaraldehyde solution is a 10% (weight/volume ratio) solution. The maximum 10% glutaraldehyde solution target weights for each size are: 0.63 grams (2mL syringe), 1.10 grams (5mL syringe), and 2.16 grams (10mL syringe). Based on these targets, the maximum quantity of glutaraldehyde coming into contact with the patient when using a single device is 0.06 grams (2mL syringe), 0.11 grams (5mL syringe), and 0.22 grams (10mL syringe) for each configuration.

## INDICATIONS / INTENDED PURPOSE

BioGlue Surgical Adhesive is indicated for use as an adjunct to standard methods of surgical repair (such as sutures, staples, and/or patches) to adhere, seal, and/or reinforce soft tissue. Indicated soft tissues are cardiac, vascular, pulmonary, and dural.

## TARGET PATIENT POPULATION

Adult patients undergoing surgery who require an adjunct to standard method of surgical repair to bond, seal, and/or reinforce cardiac, vascular, dural, and pulmonary tissue.

## INTENDED USERS

BioGlue Surgical Adhesive is intended to be used by healthcare professionals, such as surgeons, qualified in the appropriate indication.

## PERFORMANCE CHARACTERISTICS

- BioGlue reaches full bonding strength within 2 minutes.<sup>1</sup>
- BioGlue Surgical Adhesive forms strong covalent bonds with tissue and mechanically interlocks with synthetic graft material. These strong bonds lead to BioGlue having a recorded burst strength of at least 560mmHG in *in vitro* studies.<sup>1,3,4,5</sup>

### In large vascular and cardiac procedures where BioGlue was used:

- Decrease in anastomotic bleeding when compared to standard repair technique.<sup>12</sup>

### In pulmonary, large vascular and cardiac procedures where BioGlue was used:

- Reduction in ICU and hospital stays when compared to standard surgical technique.<sup>12,19,23</sup>

### In aortic dissection repair where BioGlue was used:

- Fewer pledgets, hemostatic agents, and make-up stitches were required when compared standard surgical technique.<sup>17</sup>
- Operating room time, cross clamp time, circulatory arrest time, bypass time were reduced when compared to surgical technique.<sup>18,19</sup>
- The use of platelets, plasma, and blood cells were reduced when compared to standard surgical technique.<sup>17,19</sup>

### In pulmonary procedures where BioGlue has been used:

- BioGlue has been shown to be effective in sealing air leaks when applied to a deflated or inflated lung.<sup>6,7,8,9</sup>
- Air leak duration was reduced when compared to standard surgical repair.<sup>6,8,23</sup>

### In dural procedures where BioGlue has been used:

- CSF leaks were reduced when compared to standard surgical repair.<sup>20,21,22</sup>

BioGlue Surgical Adhesive Accessories aid in the delivery of BioGlue.

## CLINICAL BENEFITS

In large vascular, cardiac, pulmonary, and dural procedures where BioGlue was used:

- Pseudoaneurysm rates have been shown to be lower when compared to published literature on standard repair technique.<sup>11,14,15</sup>
- Complication rates have been shown to be lower when compared to published literature on standard repair technique.<sup>12,2,23, 18</sup>
- Mortality rates have been shown to be lower when compared to reported literature on standard repair technique.<sup>10,11,12,13,16</sup>

## CONTRAINDICATIONS

BioGlue is contraindicated for use in cerebrovascular repairs and any intra-luminal areas.

BioGlue is not for patients with known sensitivity to materials of bovine origin.

## WARNINGS

- Do not use BioGlue as a substitute for sutures or staples in tissue approximations.
- Do not use BioGlue in a manner that would contact or obstruct circulating blood flow during or after application. BioGlue entering the circulation can result in local or embolic vascular obstruction.

- Do not use BioGlue in a manner that would obstruct circulating air or other luminal fluid flow during or after application.
- Avoid contact with nerves, eyes, or other tissue not intended for application.
- An animal study<sup>24</sup> has shown that direct application of BioGlue to the exposed phrenic nerve can cause acute nerve injury. A separate animal study<sup>25</sup> has shown that direct application of BioGlue to the surface of the sinoatrial node (SAN) of the heart can cause coagulation necrosis that extends into the myocardium, which could reach underlying conduction tissue and may cause acute, focal SAN degeneration. Subsequent animal studies<sup>26,27</sup> have demonstrated that chlorhexidine gluconate gel can protect the phrenic nerve, the myocardium, and the underlying SAN from potential injury from BioGlue use.
- Do not use BioGlue if staff is not adequately protected (e.g., wearing gloves, mask, protective clothing, and safety glasses). Unreacted glutaraldehyde may cause irritation to eye, nose, throat, or skin; induce respiratory distress; and cause local tissue necrosis. Prolonged exposure to unreacted glutaraldehyde may cause central nervous system or cardiac pathology. If contact occurs, flush affected areas immediately with water and seek medical attention.
- Polymerized BioGlue has space occupying properties. When used improperly, or applied incorrectly, serious adverse events have been reported relating to compression of adjacent anatomic structures. BioGlue should be used only when complete visualization of the target application location is possible, when it is properly primed to achieve optimal viscosity, and a minimal amount is used. Please see the *Indications/Intended Purpose and Directions for Use* sections on this label.
- Minimize use of BioGlue in patients with abnormal calcium metabolism (e.g. chronic renal failure, hyperparathyroidism). Glutaraldehyde-treated tissue has an enhanced propensity for mineralization. Laboratory experiments indicate that unreacted glutaraldehyde may have mutagenic effects.
- Do not use BioGlue in the presence of infection and use with caution in contaminated areas of the body.
- Exercise caution with repeat exposure of BioGlue in the same patient. Hypersensitivity reactions are possible upon exposure to BioGlue. Sensitization has been observed in animals.
- BioGlue contains a material of animal origin, which may be capable of transmitting infectious agents.
- The use of BioGlue in pregnant/breast feeding women has not been studied.
- BioGlue syringe and BioGlue accessories are single-use devices and should not be used in more than one patient.
- The use of BioGlue in pediatric patients has not been studied. BioGlue should not be applied circumferentially to tissue and may not allow that tissue to grow or expand.
- The suitability of BioGlue for treatment of Bronchopleural Fistula (BPF) or Lymphatic Leakage is not confirmed by sufficient data.
- When BioGlue is used in conjunction with any other material, the instructions for both products should be carefully reviewed and adhered to.

## PRECAUTIONS

- It is recommended that surgical gloves, sterile gauze pads/towels, and surgical instruments be maintained moist to minimize the potential for BioGlue inadvertently adhering to these surfaces.

- BioGlue syringe, applicator tips, spreader tips, and syringe extender tips are for single patient use only. Do not re-sterilize.
- Do not use if packages have been opened or damaged.
- Take care not to spill contents of the syringe.
- Do not compress the syringe plunger while attaching it to the syringe.
- Do not apply BioGlue in a surgical field that is too wet. This may result in poor adherence.
- Avoid tissue contact with material expelled from applicator during priming.
- BioGlue polymerizes rapidly. Priming must occur quickly, followed immediately by the application of BioGlue. Pausing between priming and application can cause polymerization within the applicator tip.
- Do not use blood saving devices when suctioning excess BioGlue from the surgical field.
- Clamp and depressurize vessels prior to applying BioGlue to targeted anastomoses.
- To prevent the entrance of BioGlue into the cardiovascular system, avoid any negative pressure during application and polymerization of BioGlue. For example, left ventricular vents should be turned off prior to the application of BioGlue. There have been reports of BioGlue being suctioned into the aorta and impeding heart valve function when used in conjunction with an active left ventricular vent.
- Do not peel away BioGlue from an unintended site, as this could result in tissue damage.
- Do not implant BioGlue into closed anatomic locations that are in immediate proximity to nerve structures.
- Due to clinical reports<sup>28</sup> of ineffective sealing when BioGlue is used in the translabrynthine approach for acoustic neuroma repairs, its use with this surgical approach is not recommended. Successful use of the product using the middle fossa or retrosigmoid approach for acoustic neuroma repair has been described in the literature and is recommended.<sup>26</sup>
- Published human clinical data has shown that excessive application of BioGlue in lung surgery can cause residual air space and atelectasis.<sup>7</sup>
- This device should be handled and disposed of in accordance with all applicable regulations including, without limitation, those pertaining to human health and safety and the environment.

#### UNDESIRABLE SIDE EFFECTS / ADVERSE INCIDENTS - OBSERVED AND POTENTIAL

Knowledge of all possible complications of soft tissue repair surgery is necessary for physicians performing these procedures. Complications specific to these types of surgeries may occur at any time during or after the procedure.

#### Observed Undesirable Side Effects / Adverse Events:

Adverse events observed during the clinical studies included the following: BioGlue applied to nontargeted tissue, Failure of BioGlue to adhere, Death, Vessel rupture and hemorrhage, Cerebrospinal fluid leak, Infection, Inflammatory, Immune systemic allergic reaction, Irreversible morbidity, Ischemia, Myocardial infarction, Neurological deficit, Organ system failure, Paraplegia, Pleural Effusion, Renal Dysfunction/failure, Respiratory Dysfunction/Failure, Stroke or Cerebral Infarction, Thromboembolism, and Thrombosis.

#### Potential Undesirable Side Effects / Adverse Events that May Occur from the Use of BioGlue:

Complications specific to the adjunctive use of BioGlue Surgical Adhesive during soft tissue repair surgery may include, but are not limited to, the following: hypersensitivity reaction such as swelling or edema at the application site, failure of product to adhere to tissue, application of adhesive to tissue not

targeted for procedure, inflammatory and immune response, allergic reaction, mineralization of tissue, local tissue necrosis, vessel obstruction, bronchial or luminal obstruction, thrombosis and thromboembolism, pulmonary emboli, injury to normal vessels or tissue, stenosis, seroma, pseudoaneurysm, and possible transmission of infectious agents from material of animal origin.

#### PACKAGING AND STORAGE

The BioGlue syringe and applicator tips are supplied sterile. Discard any unused material from opened or damaged product.

The BioGlue solutions are contained within a capped, double-chambered sterile syringe. Polymerized BioGlue is non-pyrogenic. Store below 25°C, but do not freeze.

#### ENVIRONMENTAL CONDITIONS

BioGlue Surgical Adhesive is MR safe (i.e., an item that poses no hazards in all MR environments).

#### DIRECTIONS FOR USE

Apply BioGlue Surgical Adhesive prophylactically or after a leak is detected.

#### Device Preparation

The BioGlue Surgical Adhesive Syringe delivery system consists of: syringe, syringe plunger and applicator tip.

Inside the BioGlue Syringe box there are two separate pouches. One contains the syringe and the syringe plunger, and one contains four applicator tips.

The 10mL BioGlue Syringe box includes an additional pouch containing three 12mm Spreader Tips. Visually inspect all pouches prior to use. If any breaches in the sterile barrier system are present, do not use.

1. Remove the syringe, syringe plunger, and applicator tips from their packaging. While holding the syringe upright, tap the syringe until the air bubbles in the solutions rise to the top of the syringe.


Figure 1


Figure 2

**NOTE:** Continue to hold the syringe upright during the entire assembly of the delivery system to keep the bubbles toward the top of the syringe.

2. Remove an applicator tip from its packaging and inspect the collar portion of the tip to ensure that the pointer portion is directly over the larger port. If not, rotate the locking collar on the shaft until the pointer is over the larger port. While firmly grasping the syringe, nose upward, turn the cap 90° counterclockwise and remove the cap by rocking it from side-to-side. Align the tip with the syringe using the corresponding notches on each and place the tip on the syringe.

Figure 3

**CAUTION:** ⚠ Take care not to spill solution from the syringe during assembly.

3. Lock the applicator tip in place by pushing the tip firmly toward the syringe and rotating the tip collar 90° clockwise.

Figure 4

While keeping the syringe upright, align the small and large barrels of the syringe with the corresponding syringe plunger heads and slide the plunger into the back of the syringe until resistance is felt. The syringe delivery device is now assembled.


Figure 5

**CAUTION:** ⚠ Do not lay the assembled device on its side until all air has been purged (see next paragraph).

**CAUTION:** ⚠ Before using BioGlue in the procedure, the syringe must be purged of the residual air space and the applicator tip must be primed. Refer to [Site Preparation, Syringe Air Space Removal and Applicator Tip Priming](#).

4. If using an applicator tip with a flexible extension, a desired angle may be created by bending the extension at the appropriate location to the desired angle and holding for 3-5 seconds. The angle created should be maintained for up to 5 minutes.
5. To remove occluded applicator tips, grasp the applicator tip collar, rotate the tip collar counterclockwise, and lift the tip off the syringe by rocking it side to side.

#### Site Preparation, Syringe Air Space Removal and Applicator Tip Priming

1. The target surgical field must be properly prepared prior to either removal of the residual air space, priming, or applying BioGlue. BioGlue works best when the target surgical field is dry. A dry surgical field can be described as a field that does not re-stain with blood within 4-5 seconds after wiping dry with a surgical sponge.

**CAUTION:** ⚠ Do not attempt to apply BioGlue to a field that is too wet. Application of BioGlue into a wet field may result in the failure of BioGlue to adhere.

2. The residual syringe air space must be removed prior to BioGlue application. Again, it is important to hold the assembled syringe upright to ensure that the air bubbles in the solutions are located at the top of the syringe. Purging of the air space can now be accomplished using two different methods:
  - a. Compress the plunger only until the solutions are even with the top of the syringe body. Once the residual air space has been removed the syringe is ready for priming (refer to Step 3) and immediate use.
  - b. Compress the plunger until both solutions can be visibly seen in the base of the tip. The airspace has now been removed, but this tip is now occluded with polymerized

BioGlue and will need to be changed prior to priming (refer to Step 3) and application to the target site.


Figure 6

**NOTE:** Each syringe only needs to be purged of residual airspace upon its initial use.

- Each applicator tip must be primed prior to BioGlue application. Priming ensures the BioGlue solutions are properly mixed. The surgeon should compress the plunger and expel a narrow ribbon of BioGlue approximately 3 cm long onto a sterile disposable surface (e.g., sponge, gauze, or towel).
- The surgeon should examine the material expelled during priming and ensure that it is of uniform light yellow to amber color and that it is free from air bubbles. If this material looks colorless or contains bubbles, repeat the prime as outlined in Step 2 until the device delivers a uniform liquid with no bubbles.

**CAUTION:** ⚠️ Avoid direct contact with material expelled during priming.

**CAUTION:** ⚠️ If there is evidence of syringe breakage or leakage, discard the device and open/use a new one.

- When the applicator tip has been properly primed, proceed immediately to application.

**CAUTION:** ⚠️ BioGlue polymerizes very quickly. The surgeon must apply BioGlue immediately after priming.

Pausing between priming and application can cause polymerization of BioGlue within the applicator tip. Should this occur, replace the obstructed tip with a new tip and repeat the steps for applicator tip priming. Do not continue to apply pressure to the plunger once the tip has occluded.

#### General Techniques for the Use of BioGlue in Surgery<sup>6,23,29,30</sup>

Before using BioGlue Surgical Adhesive, surgeons should become familiar by appropriate training with the surgical techniques and variations of their specific procedures. The use of BioGlue Surgical Adhesive should be practiced with the product prior to initial use in the surgical suite.

- The patient should be prepared and draped according to the hospital's standard procedures. Procedures such as entry of the chest or pleural space, cardiopulmonary bypass, clamping, and myocardial protection should follow the surgeon's standard techniques.
- The tissue surrounding the surgical site can be protected from the undesired application of BioGlue Surgical Adhesive by placing moist sterile gauze pads in these areas. Directly after application, remove gauze while adhesive is still soft, wiping away excess adhesive residue from around the site.

**CAUTION:** ⚠️ Do not use blood saving devices when suctioning excess BioGlue.

**CAUTION:** ⚠️ Clamp and depressurize vessels prior to applying BioGlue to targeted anastomoses.

**CAUTION:** ⚠️ Avoid suctioning BioGlue into the vessels when applying it to targeted anastomoses.

- If BioGlue Surgical Adhesive does adhere to an undesired location, allow the adhesive to polymerize and then gently dissect the adhesive away from the unintended area with forceps and scissors. Do not attempt to peel away the BioGlue Surgical Adhesive as this could result in tissue damage at the application site.
- For vessel repair, apply an even adhesive coating 1.2 - 3.0 mm thick for anastomosis of vessels/grfts greater than 2.5 cm in diameter; apply an even adhesive coating 0.5 - 1.0 mm for vessels/grfts less than 2.5 cm in diameter.
- The area of adhesive application should NOT be compressed or subjected to extra pressure. BioGlue works optimally when it is allowed to polymerize without any manipulation for a full two minutes. Once the adhesive has polymerized, secure with sutures as necessary.
- After adhesive polymerizes, trim away excess or irregular adhesive edges with scissors and irrigators.

#### Specific Techniques for the Use of BioGlue in Aortic Dissection Surgery<sup>6,11,13-19</sup>

- The dissected layers of the aorta should be initially cleared of blood and thrombus material and should be dried, to the extent possible, with surgical sponges.
- For the distal end of the dissection repair, insert a balloon catheter into the true lumen to define the distal terminus for the application of BioGlue. In addition, the dissected layers of the aorta should be closely approximated by inserting a dilator, sponge, or catheter into the true lumen to preserve the natural architecture of the vessel.

BioGlue should then be dispensed into the false lumen as far distally as the distal balloon catheter will allow. Filling the false lumen should proceed from distal to proximal with a spiraling out motion for smooth application. Completely fill the false lumen with BioGlue; avoid overfilling the false lumen and spilling BioGlue into the true lumen or surrounding tissue.

- For the proximal end of the dissection repair, the dissected layers of the aorta should also be closely approximated by using a dilator, sponge, or catheter. If necessary, moist gauze pads should be placed over the aortic valve leaflets to protect them from inadvertent application of BioGlue. BioGlue should then be dispensed to fill the false lumen.

Graft material may be sutured directly onto the tissues adhered and reinforced with BioGlue at both the proximal and distal aspects of the dissection repair. Allow BioGlue to completely polymerize without any manipulation for a full two minutes prior to suturing through the adhered tissue layers.

**CAUTION:** ⚠️ In order to preserve the patency of the coronary lumen in the event of dissection extension, placement of a catheter into the coronary ostia prior to BioGlue application should be considered.

#### Use of BioGlue in Lung Surgery<sup>6,9</sup>

BioGlue has been shown to be effective when applied to a deflated or inflated lung.

#### Disposal Instructions

Discard any unused material from open or damaged product by placing it in a biohazard bin.

#### Summary of Safety and Clinical Performance

<https://ec.europa.eu/tools/eudamed>  
BUDI-DI: 87723400BG35007W

#### Lifetime of BioGlue

BioGlue Surgical Adhesive is intended for long-term use (for more than 30 days). BioGlue degrades via proteolysis; it can be slow to resorb dependent on the quantity of adhesive applied and vascularity of the target tissue.

#### Reporting of Serious Incidents

Serious Incidents that occur in relation to BioGlue should be reported to the manufacturer and the competent authority of the Member State in which the patient is established. Contact information for reporting to manufacturer is below:

**Artivion, Inc.** 1655 Roberts Blvd, NW  
Kennesaw, Georgia 30144 • USA  
Phone: +1 (770) 419-3355  
Fax: +1 (770) 590-3753  
E-mail: [fieldassurance@artivion.com](mailto:fieldassurance@artivion.com)

#### REFERENCES

References regarding the information in this insert are available upon request.

<sup>1</sup>Artivion Data on File. Val-00097: BioGlue Manufacturing Process Validation

<sup>2</sup>Miscusi M, Polli FM, Forcato S, Coman MA, Ricciardi L, Ramieri A, Racco A. The use of surgical sealants in the repair of dural tears during non-instrumented spinal surgery. Eur Spine J. 2014 Aug;23(8):1761-6.

<sup>3</sup>Glock Y, et al. Experimental technique of aortoprosthetic anastomoses by gluing (BioGlue® Artivion). Presented at the Laparoscopic Aortic Aortic Surgery for Occlusive Disease and Aneurysms in Marseilles, France. January 28, 2000

<sup>4</sup>Murdoch M, et al. Cytocompatibility and mechanical properties of surgical sealants for cardiovascular applications. J Thorac Cardiovasc Surg 2018;1-8

<sup>5</sup>Gundry SR, Black K, Izutani H. Sutureless coronary artery bypass with biologic glued anastomoses: preliminary in vivo and in vitro results. J Thorac Cardiovasc Surg. 2000 Sep;120(3):473-7. doi: 10.1067/mtc.2000.108596. PubMed PMID: 10962407.

<sup>6</sup>Potaris K, et al. Experience with an albumin glutaraldehyde tissue adhesive in sealing air leaks after bullectomy. Heart Surgery Forum 2003;6(5):429-433

<sup>7</sup>Potaris K, et al. Preliminary results with the use of an albumin-glutaraldehyde tissue adhesive in lung surgery. Medical Science Monitor 2003;9(7):179-83

<sup>8</sup>Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. J Thorac Cardiovasc Surg 2006;132:105-12

<sup>9</sup>Herget GW, et al. Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. European J Cardio-Thorac Surg 2001;19(1):4-9.

<sup>10</sup>Passage J, et al. BioGlue Surgical Adhesive – An appraisal of its Indications in Cardiac Surgery. Ann Thorac Surg 2002; 74:432-7

<sup>11</sup>Fehrenbacher J, et al. Use of BioGlue in Aortic Surgery: Proper Application Techniques and results in 92 patients. The Heart Surgery Forum 2006; 9(5)

<sup>12</sup>Coselli J, et al. Prospective Randomized Study of a Protein-Based Tissue Adhesive Used as a Hemostatic and Structural Adjunct in Cardiac and Vascular Anastomotic Repair Procedures. J Am Coll Surg 2003;197:243-253

<sup>13</sup>Bavaria J, et al. Advances in the Treatment of Acute Type A Dissection: An Integrated Approach. Ann Thorac Surg 2002; 74:S1848-52

<sup>14</sup>Ma W, et al. Does BioGlue contribute to anastomotic pseudoaneurysm after thoracic aortic surgery? J Thorac Dis 2017

<sup>15</sup>Weiner J, et al. Role of Bovine Serum Albumin-Glutaraldehyde Glue in the Formation of Anastomotic Pseudoaneurysms. J Card Surg 2011; 26:76-81

<sup>16</sup>Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. Ann Thorac Surg 2002;73:707-13.

<sup>17</sup>Internal Clinical Study Report – BioGlue Surgical Adhesive Effectiveness and Safety Trial as a Surgical Adjunct in the Repair of Type A Aortic Dissection. Protocol BG1001

<sup>18</sup>Feier H, et al. The Influence of Albumin/glutaraldehyde Sealant in Early Results After Acute Type A Aortic Dissection. REV.CHIM.(Bucharest) - 70 - No. 6 – 2019

<sup>19</sup>Chao HH, Torchiana DF. BioGlue: Albumin/ Glutaraldehyde Sealant in Cardiac Surgery. J Card Surg. 2003;18:500-3

<sup>20</sup>Dusick JR, Mattozo CA, Esposito F, Kelly DF. BioGlue for prevention of postoperative cerebrospinal fluid leaks in transphenoidal surgery: A case series. Surg Neurol. 2006 Oct;66(4):371-6; discussion 376

<sup>21</sup>Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transphenoidal procedures. J Clin Neurosurg 2003; 10:92-95

<sup>22</sup>Kumar A, Maartens NF, Kaye AH. Evaluation of the use of BioGlue in neurosurgical procedures. J Clin Neurosci. 2003 Nov;10(6):661-

<sup>23</sup>Allama AM, Abd Elaziz ME. Using tissue glues for decreasing alveolar air leak in thoracic surgery. Asian Cardiovasc Thorac Ann. 2019 Apr 2.

<sup>24</sup>LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. Presented at the 34th Association of Academic Surgery 2000

<sup>25</sup>LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. Presented at Symposium for Cardiac and General Thoracic Surgery Teams. (Abstract) 2000

<sup>26</sup>LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001

<sup>27</sup>LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001

<sup>28</sup>Sen A, et al. Cerebrospinal fluid leak rate after the use of BioGlue in translabyrinthine vestibular schwannoma surgery: A prospective study. Otol Neurotol 2006;27:102-5.

<sup>29</sup>Raanani E, et al. Use of "BioGlue" in aortic surgical repair. Ann Thorac Surg 2001;72:638-40

<sup>30</sup>Coselli JS, et al. Thoracic aortic anastomoses. Operative Techniques in Thoracic and Cardiovascular Surg 2000;5(4):259-76

#### **DISCLAIMER OF WARRANTIES; LIMITS OF LIABILITY**

**ARTIVION DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES WITH RESPECT TO THIS SURGICAL ADHESIVE, INCLUDING BUT NOT LIMITED TO THE EXPRESS AND IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL ARTIVION BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES.** In the event that such disclaimer is found invalid or unenforceable for any reason: (i) any action for breach of warranty must be commenced within one year after any such claim or cause of action accrued and (ii) the remedy for any such breach is limited to the replacement of the product.

*Artivion and BioGlue are trademarks of Artivion, Inc.*

©2023 Artivion, Inc. All rights reserved. Printed in the USA.

#### **PRODUCT INFORMATION DISCLOSURE**

Handling and storage of this device by the user, as well as factors related to the patient, the diagnosis, treatment, surgical procedures, and other matters beyond manufacturer's control, may directly affect this device and the results obtained from its use.

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<sup>2</sup>Miscusi M, Polli FM, Forcato S, Coman MA, Ricciardi L, Ramieri A, Raco A. The use of surgical sealants in the repair of dural tears during non-instrumented spinal surgery. Eur Spine J. 2014 Aug;23(8):1761-6.

<sup>3</sup>Glock Y, et al. Experimental technique of aortoprotective anastomoses by gluing (BioGlue® Artivion). Presented at the Laparoscopic Aortic Surgery for Occlusive Disease and Aneurysms in Marseilles, France, January 28, 2000

<sup>4</sup>Murdoch M, et al. Cyto compatibility and mechanical properties of surgical sealants for cardiovascular applications. J Thorac Cardiovasc Surg 2018:1-8

<sup>5</sup>Gundry SR, Black K, Izutani H. Sutureless coronary artery bypass with biologic glued anastomoses: preliminary in vivo and in vitro results. J Thorac Cardiovasc Surg. 2000 Sep;120(3):473-7. doi: 10.1067/jmtc.2000.108596. PubMed PMID: 10962407.

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<sup>9</sup>Hergert GW, et al. Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. European J Cardio-Thorac Surg 2001;19(1):4-9.

<sup>10</sup>Passage J, et al. BioGlue Surgical Adhesive – An appraisal of Its Indications in Cardiac Surgery. Ann Thorac Surg 2002; 74:432-7

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<sup>12</sup>Coselli J, et al. Prospective Randomized Study of a Protein-Based Tissue Adhesive Used as a Hemostatic and Structural Adjuvant in Cardiac and Vascular Anatomic Repair Procedures. J Am Coll Surg 2003;197:243-253

- <sup>13</sup>Bavaria J, et al. Advances in the Treatment of Acute Type A Dissection: An Integrated Approach. *Ann Thorac Surg* 2002; 74:S1848-52
- <sup>14</sup>Ma W, et al. Does BioGlue contribute to anastomotic pseudoaneurysms after thoracic aortic surgery? *J Thorac Dis* 2017
- <sup>15</sup>Weiner J, et al. Role of Bovine Serum Albumin-Glutaraldehyde Glue in the Formation of Anastomotic Pseudoaneurysms. *J Card Surg* 2011; 26:76-81
- <sup>16</sup>Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. *Ann Thorac Surg* 2002;73:707-13.
- <sup>17</sup>Internal Clinical Study Report – BioGlue Surgical Adhesive Effectiveness and Safety Trial as a Surgical Adjunct in the Repair of Type A Aortic Dissection. Protocol BG1001
- <sup>18</sup>Feier H, et al. The Influence of Albumin/glutaraldehyde Sealant in Early Results After Acute Type A Aortic Dissection. *REV.CHIM.(Bucharest) - 70 - No. 6 – 2019*
- <sup>19</sup>Chao HH, Torchiana DF. BioGlue: Albumin/Glutaraldehyde Sealant in Cardiac Surgery. *J Card Surg.* 2003;18:500-3
- <sup>20</sup>Dusick JR, Mattozo CA, Esposito F, Kelly DF. BioGlue for prevention of postoperative cerebrospinal fluid leaks in transsphenoidal surgery: A case series. *Surg Neurol.* 2006 Oct;66(4):371-6; discussion 376
- <sup>21</sup>Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transsphenoidal procedures. *J Clin Neurosurg* 2003; 10:92-95
- <sup>22</sup>Kumar A, Maartens NF, Kaye AH. Evaluation of the use of BioGlue in neurosurgical procedures. *J Clin Neurosci.* 2003 Nov;10(6):661-4
- <sup>23</sup>Allama AM, Abd Elaziz ME. Using tissue glues for decreasing alveolar air leak in thoracic surgery. *Asian Cardiovasc Thorac Ann.* 2019 Apr 2.
- <sup>24</sup>LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. Presented at the 34th Association of Academic Surgery 2000
- <sup>25</sup>LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. Presented at Symposium for Cardiac and General Thoracic Surgery Teams. (Abstract) 2000
- <sup>26</sup>LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001
- <sup>27</sup>LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001
- <sup>28</sup>Sen A, et al. Cerebrospinal fluid leak rate after the use of BioGlue in translabrynthine vestibular schwannoma surgery: A prospective study. *Otol Neurotol* 2006;27:102-5.
- <sup>29</sup>Raanani E, et al. Use of "BioGlue" in aortic surgical repair. *Ann Thorac Surg* 2001;72:638-40
- <sup>30</sup>Coselli JS, et al. Thoracic aortic anastomoses. *Operative Techniques in Thoracic and Cardiovascular Surg* 2000;5(4):259-76

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<sup>26</sup>LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. [O Gel de Gluconato de Clorexidina protege os nervos expostos durante a aplicação do Adesivo Cirúrgico BioGlue.] Presented at the 33rd Association of Academic Surgery Meeting. [Apresentado no 35º Encontro da Association of Academic Surgery.] (Abstract) 2001

<sup>27</sup>LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. [O Gel de Gluconato de Clorexidina protege o miocárdio e o nó sinotrial durante a aplicação do Adesivo Cirúrgico BioGlue.] Presented at the 35th Association of Academic Surgery Meeting. [Apresentado no 35º Encontro da Association of Academic Surgery.] (Abstract) 2001

<sup>28</sup>Sen A, et al. Cerebrospinal fluid leak rate after the use of BioGlue in translabrynthine vestibular schwannoma surgery: A prospective study. [Taxa de vazamento do líquido cefalorraquidiano após o uso de BioGlue na cirurgia de schwannoma vestibular translabirintina: Um estudo prospectivo.] *Otol Neurotol* 2006;27:102-5.

<sup>29</sup>Raanani E, et al. Use of "BioGlue" in aortic surgical repair. [Uso de "BioGlue" no reparo cirúrgico da aorta.] *Ann Thorac Surg* 2001;72:638-40

<sup>30</sup>Coselli JS, et al. Thoracic aortic anastomoses. [Anastomoses da aorta torácica.] *Operative Techniques in Thoracic and Cardiovascular Surg* 2000;5(4):259-76 [Técnicas Operatórias em Cirurgia Torácica e Cardiovascular]

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*Artivion ve BioGlue, Artivion, Inc.'ticari markalarıdır.*
*©2023 Artivion, Inc. Tüm hakları saklıdır. ABD'de basılmıştır.*