(en) ENGLISH

Instructions for Use

ARTIVION™ | BioGlue®

Accessories

www.artivion.com/eifu/bioglue

A printed version of the Instructions for Use will be provided within seven days upon request to BioGlue customer service through any of the contact means listed below.

Phone: 888.427.9654 • Fax: 770.590.3753

Email: customer_service.us@artivion.com

ARTIVION™

1655 Roberts Boulevard, NW
Kennesaw, Georgia 30144, USA

Phone: +1 (888) 427-9654

FAX +1 (770) 590-3753

www.artivion.com

JOTEC GmbH

Lotzenäcker 23


D-72379 Hechingen, Germany

CE 0124

L09419.000 (2022-11-01)

EXPLANATION OF SYMBOLS

 <p>Manufacturer</p>	 <p>Temperature limit</p>
 <p>Date of manufacture</p>	 <p>Consult instructions for use / Consult electronic instructions for use</p>
 <p>EU Authorized Representative</p>	 <p>Caution</p>
 <p>Importer</p>	<p>RxONLY</p> <p>Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician</p>
 <p>Medical Device</p>	 <p>Non-pyrogenic</p>
 <p>Do not use if package is damaged</p>	 <p>Catalogue number</p>
 <p>Double sterile barrier system</p>	 <p>Batch code</p>
 <p>Sterilized using ethylene oxide</p>	 <p>Unique device identification</p>
 <p>Sterilized using irradiation</p>	 <p>Use-by date</p>
 <p>Do not re-sterilize</p>	 <p>Country of manufacture</p>
 <p>Do not re-use</p>	 <p>Quantity</p>

The following BioGlue Surgical Adhesive accessories are sold separately to aid in the delivery of BioGlue Surgical Adhesive:

Product Code	Product Description
BGAT-SY	Applicator tip
BGAT-10-SY	Syringe extender tip – 10cm
BGAT-27-SY	Syringe extender tip – 27cm
BGST-12	Spreader applicator tip – 12mm
BGST-16	Spreader applicator tip – 16mm
BGDTE-10	Delivery tip extension – 10cm
BGDTE-27	Delivery tip extension – 27cm
BGDTE-35	Delivery tip extension – 35cm

INTENDED USE

BioGlue accessories (applicator tips, spreader applicator tips, and syringe extender tips) are intended for use to deliver BioGlue Surgical Adhesive (hereinafter BioGlue) to the surgical site. Please refer to the BioGlue IFU for additional information such as indications for use, contraindications, warnings, and precautions.

PACKAGING AND STORAGE

The accessory tips are double pouched (inner pouch: Tyvek®/Tyvek® and outer pouch: Tyvek®/Poly) with 4 tips per pouch (Spreader applicator tip pouches only contain 3 tips per pouch) and 10 pouches per box. These tips are supplied sterile for single patient use only. Discard damaged or unused product components from open or compromised packaging. Product should be stored below 25°C, but do not freeze.

DIRECTIONS FOR USE

BioGlue Surgical Adhesive Preparation.

Please refer to and follow the detailed instructions in the BioGlue Surgical Adhesive IFU for Device Preparation, Site Preparation, and Syringe Air Space Removal.

Figure 1

NOTE: Continue to hold the syringe upright during the entire assembly of the BioGlue delivery system and syringe air space removal.

Attaching the Applicator Tip, Syringe Extender Tip, or Spreader Applicator Tip

1. Remove a tip from its packaging and inspect the collar portion of the tip to ensure that the pointer portion is directly over the larger port. If not, rotate the locking collar on the shaft until the pointer is over the larger port. While firmly grasping the syringe, nose upward, turn the cap 90° counterclockwise and remove the cap by rocking it from side-to-side. Align the tip with the syringe using the corresponding notches on each and place the tip on the syringe.


Figure 2


Figure 3

CAUTION: ⚠ Visually inspect all pouches prior to use. If any breaches in the sterile barrier system are present, do not use.

CAUTION: ⚠ Take care not to spill solution from the syringe during assembly.

2. Lock the tip in place by pushing the tip firmly toward the syringe and rotating the tip collar 90° clockwise.

Figure 4

BioGlue Applicator Tip, Syringe Extender Tip, or Spreader Applicator Tip Preparation and Priming

1. Each accessory tip must be primed prior to BioGlue application. Priming ensures the BioGlue solutions are properly mixed.
2. If using a tip with a flexible extension, a desired angle may be created by bending the extension at the appropriate location to the desired angle and holding for 3-5 seconds. The angle created should be maintained for up to 5 minutes.
3. To remove occluded tips, grasp the tip collar, rotate the tip collar counterclockwise, and lift the tip off the syringe by rocking it side to side.

CAUTION: ⚠ Avoid direct contact with material expelled during priming.

4. When the accessory tip assembly has been properly primed, proceed immediately to application.

CAUTION: ⚠ BioGlue polymerizes very quickly. The surgeon must apply BioGlue immediately after priming. Pausing between priming and application can cause polymerization of BioGlue within the accessory tip assembly. Should this occur, replace the obstructed applicator tip assembly with a new applicator tip then repeat the steps for priming. Do not continue to apply pressure to the plunger once the tip has occluded.

Disposal Instructions

Discard any unused material from open or damaged product by placing it in a biohazard bin

WARNINGS ⚠

- BioGlue applicator tips, spreader applicator tips, and syringe extender tips are for use exclusively with BioGlue Surgical Adhesive syringes.

- BioGlue applicator tips, spreader applicator tips, and syringe extender tips are for single patient use only. Do not re-sterilize. Re-use of BioGlue Surgical Adhesive accessories can result in infection.
- This device should be handled and disposed of in accordance with all applicable regulations including, without limitation, those pertaining to human health and safety and the environment.
- Do not use if packages have been opened unintentionally or damaged.

Artivion and BioGlue are registered trademarks of Artivion, Inc. Tyvek is a registered trademark owned by affiliates of DuPont de Nemours, Inc.

©2022 Artivion, Inc. All rights reserved. Printed in the USA.

NOTE: BioGlue Surgical Adhesive accessories are not intended to be used as surgical tools (i.e. pressing, lifting, and/or holding tissues).

Observed Undesirable Side Effects / Adverse Events:

There were no observed undesirable side effects associated with BioGlue Surgical Adhesive accessories during the clinical studies.

Summary of Safety and Clinical Performance

<https://ec.europa.eu/tools/eudamed>

BUDI-DI: 87723400BG35007W

PERFORMANCE CHARACTERISTICS

BioGlue Surgical Adhesive Accessories aid in the delivery of BioGlue.

Reporting of Serious Incidents

Serious Incidents that occur in relation to BioGlue Surgical Adhesive accessories should be reported to the manufacturer and the competent authority of the Member State in which the patient is established. Contact information for reporting to manufacturer is below:

Artivion, Inc. 1655 Roberts Blvd, NW
 Kennesaw, Georgia 30144 • USA
 Phone: +1 (888) 427-9654
 Fax: +1 (770) 590-3753
E-mail: fieldassurance@artivion.com

PRODUCT INFORMATION DISCLOSURE

Handling and storage of this device by the user, as well as factors related to the patient, the diagnosis, treatment, surgical procedures, and other matters beyond manufacturer's control, may directly or indirectly affect this device and the results obtained from its use.

DISCLAIMER OF WARRANTIES; LIMITS OF LIABILITY

ARTIVION DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES WITH RESPECT TO THE APPLICATOR TIPS, THE SYRINGE EXTENDER TIPS AND THE SPREADER APPLICATOR TIPS, INCLUDING BUT NOT LIMITED TO THE EXPRESS AND IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL ARTIVION BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES. In the event that such disclaimer is found invalid or unenforceable for any reason: (i) any action for breach of warranty must be commenced within one year after any such claim or cause of action accrued and (ii) the remedy for any such breach is limited to the replacement of the product.

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©2022 Artivion, Inc. Tüm hakları saklıdır. ABD'de basılmıştır.