(en) ENGLISH

Instructions for Use

ARTIVION™ | BioGlue®

Delivery Tip Extension

www.artivion.com/eifu/bioglue

A printed version of the Instructions for Use will be provided within seven days upon request to BioGlue customer service through any of the contact means listed below.

Phone: 888.427.9654 • Fax: 770.590.3753

Email: customer_service.us@artivion.com

ARTIVION™

1655 Roberts Boulevard, NW
Kennesaw, Georgia 30144, USA

Phone: +1 (888) 427-9654

FAX +1 (770) 590-3753

www.artivion.com

EC REP

JOTEC GmbH


Lotzenäcker 23

D-72379 Hechingen, Germany

CE 0124

L09420.000 (2022-11-03)

EXPLANATION OF SYMBOLS

	Manufacturer
	Date of manufacture
	EU Authorized Representative
	Importer
	Medical Device
	Do not use if package is damaged
	Double sterile barrier system
	Sterilized using ethylene oxide
	Do not re-sterilize
	Do not re-use
	Temperature limit
	Consult instructions for use / Consult electronic instructions for use
	Caution
RxONLY	Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician
	Non-pyrogenic
	Catalogue number
	Batch code
	Unique Device Identifier
	Use-by date
	Country of manufacture
	Quantity

DEVICE DESCRIPTION

The delivery tip extension (DTE) is intended to be used with BioGlue® Surgical Adhesive (BioGlue) syringe standard applicator tips. The DTEs, available in a variety of lengths to accommodate different surgical needs, slip over the standard applicator tips to deliver BioGlue for use as stated in the BioGlue product Instructions for Use (IFU). The DTE is made from biocompatible, USP Class VI compliant polymer resins.

INTENDED USE

The DTE is intended for use with BioGlue syringe standard applicator tips to deliver BioGlue to the surgical site. Please refer to the BioGlue IFU for additional information such as indications for use, contraindications, warnings, and precautions.

PACKAGING AND STORAGE

The DTEs are double pouched (inner and outer pouch: Tyvek®/Poly) with 2 DTEs per pouch and 10 pouches per box and are supplied sterile for single patient use only. Discard damaged or unused product components from open or compromised packaging. Product should be stored at room temperature.

DIRECTIONS FOR USE

BioGlue Preparation.

Please refer to and follow the detailed instructions in the BioGlue Syringe IFU for Device Preparation Steps 1-5, Site Preparation, and Syringe Air Space Removal Steps 1 and 2.

NOTE: Continue to hold the syringe upright during the entire assembly of the BioGlue delivery system and syringe air space removal.

Attaching the Delivery Tip Extension

1. Within the sterile field, carefully remove a DTE from the packaging.

NOTE: The DTE is sold separately.

2. Insert the flared end of the DTE over the applicator tip and push until fully seated. The DTE is fully seated when the tapered end of the standard applicator tip contacts the inner tapered portion of the DTE, as shown in Figure 1.

Figure 1 – Fully seated

Figure 2 – Not fully seated

CAUTION: ⚠ Do not compress the plunger while attaching the DTE to the applicator tip. Compressing the plunger can be avoided by: (1) Removing the plunger while attaching the DTE and reinserting the plunger after the DTE is attached, or (2) Taking care not to manually compress the plunger or compress it against any surface.

Figure 3 – Wrong

Figure 4 – Right

BioGlue Applicator Tip and Delivery Tip Extension Preparation and Priming

1. Each applicator tip/DTE assembly must be primed prior to BioGlue application. Priming ensures that the BioGlue solutions are properly mixed. The surgeon should compress the plunger and expel a narrow ribbon of BioGlue approximately 3 cm long onto a sterile disposable surface (e.g., sponge, gauze or towel).
2. The surgeon should examine the material expelled during priming and ensure that it is of uniform light yellow to amber color and that it is free from air bubbles. If this material looks colorless or contains bubbles, repeat the prime as outlined in Step 1 until the device delivers a uniform liquid with no bubbles.

CAUTION: ⚠ Avoid direct contact with material expelled during priming.

3. When the applicator tip/DTE assembly has been properly primed, proceed immediately to application.

CAUTION: ⚠ BioGlue polymerizes very quickly. The surgeon must apply BioGlue immediately after priming. Pausing between priming and application can cause polymerization of BioGlue within the applicator tip/DTE assembly. Should this occur, replace the obstructed applicator tip/DTE assembly with a new applicator tip and attach another DTE; then repeat the steps for priming. Do not continue to apply pressure to the plunger once the tip has occluded.

Disposal Instructions

Discard any unused material from open or damaged product by placing it in a biohazard bin.

WARNINGS ⚠

- BioGlue DTE are for use exclusively with BioGlue Surgical Adhesive Syringe and applicator tip.
- BioGlue DTE are for single patient use only. Do not re-sterilize. Re-use of BioGlue DTE can result in infection.
- This device should be handled and disposed of in accordance with all applicable regulations including, without limitation, those pertaining to human health and safety and the environment.
- Do not use if packages have been opened unintentionally or damaged.

NOTE: The DTE is intended to be used to facilitate delivery of BioGlue and is not intended to be used as a surgical tool (i.e. pressing, lifting, and/or holding tissues).

Observed Undesirable Side Effects / Adverse Events:

There were no observed undesirable side effects associated with BioGlue Surgical Adhesive accessories during the clinical studies.

Summary of Safety and Clinical Performance

<https://ec.europa.eu/tools/eudamed>

BUDI-DI: 87723400BG35007W

PERFORMANCE CHARACTERISTICS

BioGlue Surgical Adhesive Accessories aid in the delivery of BioGlue.

Reporting of Serious Incidents

Serious Incidents that occur in relation to BioGlue Surgical Adhesive accessories should be reported to the manufacturer and the competent authority of the Member State in which the patient is established. Contact information for reporting to manufacturer is below:

Artivion, Inc. 1655 Roberts Blvd, NW • Kennesaw, Georgia 30144 • USA
Phone: +1 (888) 427-9654 • Fax: +1 (770) 590-3753
E-mail: fieldassurance@artivion.com

PRODUCT INFORMATION DISCLOSURE

Handling and storage of this device by the user, as well as the factors related to the patient, the patient's diagnosis, treatment, surgical procedures, and other matters beyond manufacturer's control, may directly or indirectly affect this device and the results obtained from its use.

DISCLAIMER OF WARRANTIES; LIMITS OF LIABILITY

ARTIVION DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES WITH RESPECT TO THIS DELIVERY TIP EXTENSION, INCLUDING BUT NOT LIMITED TO THE EXPRESS AND IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL ARTIVION BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES. In the event that such disclaimer is found invalid or unenforceable for any reason: (i) any action for breach of warranty must be commenced within one year after any such claim or cause of action accrued and (ii) the remedy for any such breach is limited to the replacement of the product.

Artivion and BioGlue are registered trademarks of Artivion, Inc. Tyvek is a registered trademark owned by affiliates of DuPont de Nemours, Inc.

©2022 Artivion, Inc. All rights reserved. Printed in the USA.

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©2022 Artivion, Inc. Tüm hakları saklıdır. ABD'de basılmıştır.