BioGlue®

Surgical Adhesive

Instructions for Use

Syringe

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RxONLY

PRODUCT DESCRIPTION

BioGlue® Surgical Adhesive is composed of bovine serum albumin and glutaraldehyde. The product's two components are mixed within the applicator tip of the syringe during application to tissue. Polymerization of the surgical adhesive begins immediately upon application, reaching bonding strength within 2 minutes.

INDICATIONS

BioGlue Surgical Adhesive is indicated for use as an adjunct to standard methods of surgical repair (such as sutures, staples, electrocautery, and/or patches) to bond, seal, and/or reinforce soft tissue. BioGlue may also be applied alone to seal and/or reinforce damaged parenchyma when other ligation or conventional procedures are ineffective or impractical. Indicated soft tissues are cardiac, vasopulmonary, genitourinary, dural, alimentary (esophageal, gastrointestinal, and colorectal), and other abdominal (pancreatic, splenic, hepatic, and biliary). Additionally, BioGlue is used in the fixation of surgical meshes in hernia repair.

CONTRAINDICATIONS

BioGlue is contraindicated for use in cerebrovascular repairs and any intra-luminal areas. BioGlue is not for patients with known sensitivity to materials of bovine origin.

WARNINGS

- Do not use BioGlue as a substitute for sutures or staples in tissue approximations.
- Do not use BioGlue in a manner that would contact or obstruct circulating blood flow during or after application. BioGlue entering the circulation can result in local or embolic vascular obstruction.
- Do not use BioGlue in a manner that would obstruct circulating air or other luminal fluid flow during or after application.
- Avoid contact with nerves, eyes, or other tissue not intended for application.
- An animal study²² has shown that direct application of BioGlue to the exposed phrenic nerve can cause acute nerve injury. A separate animal study²³ has shown that direct application of BioGlue to the surface of the sinoatrial node (SAN) of the heart can cause coagulation necrosis that extends into the myocardium, which could reach underlying conduction tissue and may cause acute, focal SAN degeneration. Subsequent animal studies²⁴⁻²⁶ have demonstrated that chlorhexidine gluconate gel (e.g., Surgilube[®]) can protect the phrenic nerve, the myocardium, and the underlying SAN from potential injury from BioGlue if used.
- Do not use BioGlue if staff is not adequately protected (e.g., wearing gloves, mask, protective clothing, and safety glasses). Unreacted glutaraldehyde may cause irritation to eye, nose, throat, or skin; induce respiratory distress, and cause local tissue necrosis. Prolonged exposure to unreacted glutaraldehyde may cause central nervous system or cardiac pathology. If contact occurs, flush affected areas immediately with water and seek medical attention.
- Do not use BioGlue in the presence of infection and use with caution in contaminated areas of the body.
- Exercise caution with repeat exposure of BioGlue in the same patient. Hypersensitivity reactions are possible upon exposure to BioGlue. Sensitization has been observed in animals.
- BioGlue contains a material of animal origin, which may be capable of transmitting infectious agents.

- Do not use BioGlue if staff is not adequately protected (e.g., wearing gloves, mask, protective clothing, and safety glasses). Unreacted glutaraldehyde may cause irritation to eye, nose, throat, or skin; induce respiratory distress, and cause local tissue necrosis. Prolonged exposure to unreacted glutaraldehyde may cause central nervous system or cardiac pathology. If contact occurs, flush affected areas immediately with water and seek medical attention.
- Do not use BioGlue in the presence of infection and use with caution in contaminated areas of the body.
- Exercise caution with repeat exposure of BioGlue in the same patient. Hypersensitivity reactions are possible upon exposure to BioGlue. Sensitization has been observed in animals.
- BioGlue contains a material of animal origin, which may be capable of transmitting infectious agents.

PRECAUTIONS

- It is recommended that surgical gloves, sterile gauze pads/towels, and surgical instruments be maintained moist to minimize the potential for BioGlue inadvertently adhering to these surfaces.
- BioGlue syringe, applicator tips, and applicator tip extenders are for single patient use only. Do not re-sterilize.
- Do not use if packages have been opened or damaged.
- Take care not to spill contents of the syringe.
- Do not compress the syringe plunger while attaching it to the syringe.
- Do not apply BioGlue in a surgical field that is too wet. This may result in poor adherence.
- Avoid tissue contact with material expelled from applicator during priming.
- BioGlue polymerizes rapidly. Priming must occur quickly, followed immediately by the application of BioGlue. Pausing between priming and application can cause polymerization within the applicator tip.
- Do not use blood saving devices when suctioning excess BioGlue from the surgical field.
- Clamp and depressurize vessels prior to applying BioGlue to targeted anastomoses.
- To prevent the entrance of BioGlue into the cardiovascular system, avoid any negative pressure during application and polymerization of BioGlue. For example, left ventricular vents should be turned off prior to the application of BioGlue. There have been reports of BioGlue being suctioned into the aorta and impeding heart valve function when used in conjunction with an active left ventricular vent.

- Do not peel away BioGlue from an unintended site, as this could result in tissue damage.
- An animal study²² has shown that the circumferential application of BioGlue restricts dilatation on growing tissue. These findings suggest caution with the circumferential use of BioGlue in children.
- Do not implant BioGlue into closed anatomic locations that are in immediate proximity to nerve structures.
- Due to clinical reports²⁰ of ineffective sealing when BioGlue is used in the [translabrynthine approach](#) for acoustic neuroma repairs, its use with this surgical approach is not recommended. Successful use of the product using the middle fossa or retrosigmoid approach for acoustic neuroma repair has been described in the literature and is recommended.⁷
- Published human clinical data has shown that excessive application of BioGlue in lung surgery can cause residual air space and atelectasis.¹⁹

ADVERSE INCIDENTS

Knowledge of all possible complications of soft tissue repair surgery is necessary for physicians performing these procedures. Complications specific to these types of surgeries may occur at any time during or after the procedure.

Complications specific to the adjunctive use of BioGlue Surgical Adhesive during soft tissue repair surgery may include, but are not limited to, the following: failure of product to adhere to tissue, application of adhesive to tissue not targeted for procedure, inflammatory and immune response, allergic reaction, mineralization of tissue, local tissue necrosis, vessel obstruction, bronchial or luminal obstruction, thrombosis and thromboembolism, pulmonary emboli, injury to normal vessels or tissue, and possible transmission of infectious agents from material of animal origin.

PACKAGING AND STORAGE

The BioGlue syringe and applicator tips are supplied sterile for single-patient use only. Discard any unused material from opened or damaged product.

The BioGlue solutions are contained within a capped, double-chambered sterile syringe. Polymerized BioGlue is non-lyophilized. Store below 25°C, but do not freeze.

DIRECTIONS FOR USE

Apply BioGlue Surgical Adhesive prophylactically or after a leak is detected.

Device Preparation

The BioGlue Surgical Adhesive Syringe delivery system consists of: syringe, syringe plunger and applicator tip.

Inside the BioGlue Syringe box there are two separate pouches. One contains the syringe and the syringe plunger and one contains four applicator tips.

The 10mL BioGlue Syringe box includes an additional pouch containing three, 12mm Spreader Tips.

- Remove the syringe, syringe plunger, and applicator tips from their packaging. While holding the syringe upright, tap the syringe until the air bubbles in the solutions rise to the top of the syringe.

Figure 1

NOTE: Continue to hold the syringe upright during the entire assembly of the delivery system to keep the bubbles toward the top of the syringe.

- Remove an applicator tip from its packaging and inspect the collar portion of the tip to ensure that the pointer portion is directly over the larger port. If not, rotate the locking collar on the shaft until the pointer is over the larger port.


Figure 2

- While firmly grasping the syringe, nose upward, turn the cap 90° counterclockwise and remove the cap by rocking it from side-to-side. Align the tip with the syringe using the corresponding notches on each and place the tip on the syringe.


Figure 3

CAUTION: Take care not to spill solution from the syringe during assembly.

- Lock the applicator tip in place by pushing the tip firmly toward the syringe and rotating the tip collar 90° clockwise.

Figure 4

- While keeping the syringe upright, align the small and large barrels of the syringe with the corresponding syringe plunger heads and slide the plunger into the back of the syringe until resistance is felt. The syringe delivery device is now assembled.

Figure 5

CAUTION: Do not lay the assembled device on its side until all air has been purged (see next paragraph).

CAUTION: Before using BioGlue in the procedure, the syringe must be purged of the residual air space and the applicator tip must be primed. Refer to [Site Preparation, Syringe Air Space Removal and Applicator Tip Priming](#).

- If using an applicator tip with a flexible extension, a desired angle may be created by bending the extension at the appropriate location to the desired angle and holding for 3-5 seconds. The angle created should be maintained for up to 5 minutes.
- To remove occluded applicator tips, grasp the applicator tip collar, rotate the tip collar counterclockwise, and lift the tip off the syringe by rocking it side to side.

Site Preparation, Syringe Air Space Removal and Applicator Tip Priming

- The target surgical field must be properly prepared prior to either removal of the residual air space, priming, or applying BioGlue. BioGlue works best when the target surgical field is dry. A dry surgical field can be described as a field that does not remain with blood within 4-5 seconds after wiping dry with a surgical sponge.

CAUTION: Do not attempt to apply BioGlue to a field that is too wet. Application of BioGlue into a wet field may result in the failure of BioGlue to adhere.

- The residual syringe air space must be removed prior to BioGlue application. Again, it is important to hold the assembled syringe upright to ensure that the air bubbles in the solutions are located at the top of the syringe. Purging of the air space can now be accomplished using two different methods:

- Compress the plunger only until the solutions are even with the top of the syringe body. Once the residual air space has been removed the syringe is ready for priming (refer to Step 3) and immediate use.
- Compress the plunger until both solutions can be visibly seen in the base of the tip. The airspace has now been removed, but this tip is now occluded with polymerized BioGlue and will

need to be changed prior to priming (refer to Step 3) and application to the target site.


Figure 6

NOTE: Each syringe only needs to be purged of residual airspace upon its initial use.

- Each applicator tip must be primed prior to BioGlue application. Priming ensures the BioGlue solutions are properly mixed. The surgeon should compress the plunger and expel a narrow ribbon of BioGlue approximately 3 cm long onto a sterile disposable surface (e.g., sponge, gauze, or towel).
- The surgeon should examine the material expelled during priming and ensure that it is of uniform light yellow to amber color and that it is free from air bubbles. If this material looks colorless or contains bubbles, repeat the prime as outlined in Step 2 until the device delivers a uniform liquid with no bubbles.

CAUTION: Avoid direct contact with material expelled during priming.

- When the applicator tip has been properly primed, proceed immediately to application.

CAUTION: BioGlue polymerizes very quickly. The surgeon must apply BioGlue immediately after priming. Pausing between priming and application can cause polymerization of BioGlue within the applicator tip. Should this occur, replace the applicator tip with a new tip and repeat the steps for priming. Do not continue to apply pressure to the plunger once the tip has occluded.

General Techniques for the Use of BioGlue in Surgery¹²

Before using BioGlue Surgical Adhesive, surgeons should become familiar by appropriate training with the surgical techniques and variations of their specific procedures. The use of BioGlue Surgical Adhesive should be practiced with the product prior to initial use in the surgical suite.

- The patient should be prepared and draped according to the hospital's standard procedures. Procedures such as entry of the chest or pleural space, cardiopulmonary bypass, clamping, and myocardial protection should follow the surgeon's standard techniques.
- The tissue surrounding the surgical site can be protected from the undesired application of BioGlue Surgical Adhesive by placing moist sterile gauze pads in these areas. Directly after application, remove gauze while adhesive is still soft, wiping away excess adhesive residue from around the site.

CAUTION: Do not use blood saving devices when suctioning excess BioGlue.

CAUTION: Clamp and depressurize vessels prior to applying BioGlue to targeted anastomoses.

CAUTION: Avoid suctioning BioGlue into the vessels when applying it to targeted anastomoses.

- If BioGlue Surgical Adhesive does adhere to an undesired location, allow the adhesive to polymerize and then gently dissect the adhesive away from the unintended area with forceps and scissors. Do not attempt to peel away the BioGlue Surgical Adhesive as this could result in tissue damage at the application site.
- For vessel repair apply an even adhesive coating 1.2 - 3.0 mm thick for anastomosis of vessels/grfts greater than 2.5 cm in diameter; apply an even adhesive coating 0.5 - 1.0 mm for vessels/grfts less than 2.5 cm in diameter.
- For parenchymal repair apply an even adhesive coating 1.5 - 3.0 mm thick.
- The area of adhesive application should NOT be compressed or subjected to extra pressure. Once the adhesive has polymerized secure with sutures as necessary.
- After adhesive polymerizes, trim away excess or irregular adhesive edges with scissors and pickups.

Specific Techniques for the Use of BioGlue in Aortic Dissection Surgery^{1,2,16}

- The dissected layers of the aorta should be initially cleared of blood and thrombus material and should be dried, to the extent possible, with surgical sponges.
- For the distal end of the dissection repair, insert a balloon catheter into the true lumen to define the distal terminus for the application of BioGlue. In addition, the dissected layers of the aorta should be closely approximated by inserting a dilator, sponge, or catheter into the true lumen to preserve the natural architecture of the vessel.

BioGlue should then be dispensed into the false lumen as far distally as the distal balloon catheter will allow. Filling the false lumen should proceed from distal to proximal with a spiraling out motion for smooth application. Completely fill the false lumen with BioGlue; avoid overfilling the false lumen and spilling BioGlue into the true lumen or surrounding tissue.

- For the proximal end of the dissection repair, the dissected layers of the aorta should also be closely approximated by using a dilator, sponge, or catheter. If necessary, moist gauze pads should be placed over the aortic valve leaflets to protect them from inadvertent application of BioGlue. BioGlue should then be dispensed to fill the false lumen.

Graft material may be sutured directly onto the tissues adhered and reinforced with BioGlue at both the proximal and distal aspects of the dissection repair. Allow BioGlue to completely polymerize without any manipulations for a full two minutes prior to suturing through the adhered tissue layers.

Use of BioGlue in Lung Surgery^{17,20}

BioGlue has been shown to be effective when applied to a deflated or inflated lung.

RxONLY	= Prescription Device
	= Latex Free
	= Non-pyrogenic
	= Do not use if package is damaged

In the United States of America Federal law restricts a prescription device to sale by or on the order of a physician.

In all other regions of the world, this device should not be used except on the order of a physician.

REFERENCES

References regarding the information in this insert are available upon request.

- Bavaria JE, et al. Advances in the treatment of acute type A dissection: an integrated approach. *Ann Thorac Surg* 2002;74:51846-52.
- Bavaria JE. Acute aortic dissection: current perspectives. *J Penn. Assoc. Thorac Surg* 1999;7-8.
- Coselli JS, et al. Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. *J Am Coll Surg* 2003;197(2):243-253.
- Coselli JS, et al. Thoracic aortic anastomoses. Operative Techniques in Thoracic and Cardiovascular Surg 2000;5(4):259-76.
- Glock et al. Experimental techniques of aorto-prosthetic anastomosis by gluing (BioGlue - CryoLife). *Laparoscopic Aortic/Ulcer Surgery for Occlusive Disease and Aneurysms*. Jan. 28, 2000.
- Hasse J, et al. Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. *European J Cardio-Thorac Surg* 2001;19(1):4-9.
- Kumar A, et al. Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (In Press) *J Clinical Neuroscience*.
- Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transphenoidal procedures. *J Clinical Neuroscience* 2003;10(1):92-95.
- Lang H, et al. Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). A pilot study. Presented at the XVIIIth European Association of Urology on 13 March 2003 (Abstract).
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. Presented at the 34th Association of Academic Surgery 2000.
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. Presented at Symposium for Cardiac and General Thoracic Surgery Teams. (Abstract) 2000.
- LeMaire SA, et al. BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. *Ann Thorac Surg* 2002;73(5):1500-5.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
- Passage J, et al. BioGlue surgical adhesive-an appraisal of its indications in cardiac surgery. *Ann Thorac Surg* 2002;74:432-7.
- Potarnik K, et al. Experience with an albumin-glutaraldehyde tissue adhesive in sealing air leaks after bullectomy. *Heart Surg Forum* 2003;6(4):429-433.
- Potarnik K, et al. Preliminary results with the use of an albumin-glutaraldehyde tissue adhesive in lung surgery. *Medical Science Monitor* 2003;9(7):179-83.
- Raanani E, et al. Use of "BioGlue" in aortic surgical repair. *Ann Thorac Surg* 2001;72:638-40.
- Sen A, et al. Cerebrospinal fluid leak rate after the use of BioGlue in transabyrinthine vestibular schwannoma surgery: A prospective study. *Otol Neurotol* 2006;27:102-5.
- Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. *J Thorac Cardiovasc Surg* 2006;132:105-12.
- Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. *Ann Thorac Surg* 2002;73:707-13.

PRODUCT INFORMATION DISCLOSURE

Handling and storage of this device by the user, as well as factors related to the patient, the diagnosis, treatment, surgical procedures, and other matters beyond manufacturer's control, may directly affect this device and the results obtained from its use. This device should not be used except on the order of a physician.

DISCLAIMER OF WARRANTIES: LIMITS OF LIABILITY

CRYOLIFE DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES WITH RESPECT TO THIS SURGICAL ADHESIVE, INCLUDING BUT NOT LIMITED TO THE EXPRESS AND IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL CRYOLIFE BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES. In the event that such disclaimer is found invalid or unenforceable for any reason: (i) any action for breach of warranty must be commenced within one year after any such claim or cause of action accrued and (ii) the remedy for any such breach is limited to the replacement of the product.

Prices, specifications, and availability are subject to change without notice.

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1. Bavaria JE, et al. Advances in the treatment of acute type A dissection: an integrated approach. *Ann Thorac Surg* 2002;74:1848-52.

2. Bavaria JE. Acute aortic dissection: current perspectives. *J Penn. Assoc. Thorac Surg* 1999;7-8.

3. Coselli JS, et al. Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. *J Am Coll Surg* 2003;197(2):243-253.

4. Coselli JS, et al. Thoracic aortic anastomoses. *Operative Techniques in Thoracic and Cardiovascular Surg* 2000;5(4):225.

5. Glock et al. Experimental techniques of aortic prosthetic anastomosis by gluing (BioGlue – CryoLife). *Laparoscopic Aortic/Thoracic Surgery for Ocular Disease and Aneurysms*. Jan. 28, 2000.

6. Hasse J, et al. Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. *European J Cardio-Thorac Surg* 2001;19(1):4-9.

7. Kumar A, et al. Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (In Press) *J Clinical Neuroscience*.

8. Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transphenoidal procedures. *J Clinical Neuroscience* 2003;10(1):92-95.

9. Lang H, et al. Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). A pilot study. Presented at the XVIIIth European Association of Urology on 13 March 2003 (Abstract).

10. LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. Presented at the 34th Association of Academic Surgery 2000.

11. LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. Presented at Symposium for Cardiac and General Thoracic Surgery Teams. (Abstract) 2000.

12. LeMaire SA, et al. BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. *Ann Thorac Surg* 2002;73(5):1500-5.

13. LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.

14. LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.

15. Passage J, et al. BioGlue surgical adhesive-an appraisal of its indications in cardiac surgery. *Ann Thorac Surg* 2002;74:432-7.

16. Potaris K, et al. Experience with an albumin-glutaraldehyde tissue adhesive in sealing air leaks after bullectomy. *Heart Surgery Forum* 2003;6(5):429-433.

17. Potaris K, et al. Preliminary results with the use of an albumin-glutaraldehyde tissue adhesive in lung surgery. *Medical Science Monitor* 2003;9(7):179-83.

18. Raanani E, et al. Use of "BioGlue" in aortic surgical repair. *Ann Thorac Surg* 2001;72:638-40.

19. Sen A, et al. Cerebrospinal fluid leak rate after the use of BioGlue in transabyrinthine vestibular schwannoma surgery: A prospective study. *Otol Neurotol* 2006;27:102-5.

20. Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. *J Thorac Cardiovasc Surg* 2006;132:105-12.

21. Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. *Ann Thorac Surg* 2002;73:707-13.

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1. Bavaria JE, et al. Advances in the treatment of acute type A dissection: an integrated approach. *Ann Thorac Surg* 2002;74:1848-52.
2. Bavaria JE, et al. Aortic dissection: current perspectives. *J Penn. Assoc. Thorac Surg* 1999;7-8.
3. Coselli JS, et al. Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. *J Am Coll Surg* 2003;197(2):243-253.
4. Coselli JS, et al. Thoracic aortic anastomoses. *Operative Techniques in Thoracic and Cardiovascular Surg* 2000;5(4):259-76.
5. Glock et al. Experimental techniques of aortic prosthetic anastomosis by gluing (BioGlue – CryoLife). *Laparoscopic Aortic/Thoracic Surgery for Occlusive Disease and Aneurysms*. Jan. 28, 2000.
6. Hasse J, et al. Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. *European J Cardio-Thorac Surg* 2001;19(1):4-9.
7. Kumar A, et al. Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (In Press) *J Clinical Neuroscience*
8. Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transphenoidal procedures. *J Clinical Neuroscience* 2003;10(1):92-95.
9. Lang H, et al. Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). A pilot study. Presented at the XVIIIth European Association of Urology on 13 March 2003 (Abstract).
10. LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. Presented at the 34th Association of Academic Surgery 2000.
11. LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. Presented at Symposium for Cardiac and General Thoracic Surgery Teams. (Abstract) 2000.
12. LeMaire SA, et al. BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. *Ann Thorac Surg* 2002;73(5):1500-5.
13. LeMaire SA, et al. Chlorhexidine Gluconate Gel increases exposed nerves during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
14. LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
15. Passage J, et al. BioGlue surgical adhesive-an appraisal of its indications in cardiac surgery. *Ann Thorac Surg* 2002;74:432-7.
16. Potaris K, et al. Experience with an albumin-glutaraldehyde tissue adhesive in sealing air leaks after bullectomy. *Heart Surgery Forum* 2003;6(5):429-433.
17. Potaris K, et al. Preliminary results with the use of an albumin-glutaraldehyde tissue adhesive in lung surgery. *Medical Science Monitor* 2003;9(7):179-83.
18. Raanani E, et al. Use of "BioGlue" in aortic surgical repair. *Ann Thorac Surg* 2001;72:638-40.
19. Sen A, et al. Cerebrospinal fluid leak rate after the use of BioGlue transabdominal vesicular schwannoma surgery: A prospective study. *Otol Neurotol* 2006;27:102-5.
20. Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. *J Thorac Cardiovasc Surg* 2006;132:105-12.
21. Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. *Ann Thorac Surg* 2002;73:707-13.

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- Bavaria JE, et al. Advances in the treatment of acute type A dissection: an integrated approach. *Ann Thorac Surg* 2002;74:S1848-52.
- Bavaria JE. Acute aortic dissection: current perspectives. *J Penn Assoc. Thorac Surg* 1999;7-8.
- Coselli JS, et al. Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. *J Am Coll Surg* 2003;197(2):243-253.
- Coselli JS, et al. Thoracic aortic anastomoses. *Operative Techniques in Thoracic and Cardiovascular Surg* 2000;5(4):259-76.
- Glock et al. Experimental techniques of aortic prosthetic anastomosis by gluing (BioGlue - CryoLife). *Laparoscopic Aortic/ail Surgery for Occlusive Disease and Aneurysms*. Jan. 28, 2000.
- Hasse J, et al. Experimental use of an albumin-glyceraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. *European J Cardio-Thorac Surg* 2001;19(1):4-9.
- Kumar A, et al. Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (In Press) *J Clinical Neuroscience*.
- Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transphenoidal procedures. *J Clinical Neuroscience* 2003;10(1):92-95.
- Lang H, et al. Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). A pilot study. Presented at the XVIIth European Association of Urology on 13 March 2003 (Abstract).
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. Presented at the 34th Association of Academic Surgery 2000.
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. Presented at Symposium for Cardiac and General Thoracic Surgery Teams. (Abstract) 2000.
- LeMaire SA, et al. BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. *Ann Thorac Surg* 2002;73(5):1500-5.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
- Passage J, et al. BioGlue surgical adhesive-an appraisal of its indications in cardiac surgery. *Ann Thorac Surg* 2002;74:432-7.
- Potaris K, et al. Experience with an albumin-glyceraldehyde tissue adhesive in sealing air leaks after bullectomy. *Heart Surgery Forum* 2003;6(5):429-433.
- Potaris K, et al. Preliminary results with the use of an albumin-glyceraldehyde tissue adhesive in lung surgery. *Medical Science Monitor* 2003;9(7):179-83.
- Raanan E, et al. Use of "BioGlue" in aortic surgical repair. *Ann Thorac Surg* 2001;72:638-40.
- Sen A, et al. Cerebrospinal fluid leak rate after the use of BioGlue in translabrynthine vestibular schwannoma surgery: A prospective study. *Otol Neurotol* 2006;27:102-5.
- Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. *J Thorac Cardiovasc Surg* 2006;132:105-12.
- Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. *Ann Thorac Surg* 2003;73:707-13.

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1. Bavaria JE, et al. Advances in the treatment of acute type A dissection: an integrated approach. *Ann Thorac Surg* 2002;74:1848-52.
2. Bavaria JE. Acute aortic dissection: current perspectives. *J Penn. Assoc. Thorac Surg* 1999;7:8.
3. Coselli JS, et al. Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. *J Am Coll Surg* 2003;197(2):243-253.
4. Coselli JS, et al. Thoracic aortic anastomoses. Operative Techniques in Thoracic and Cardiovascular Surg 2000;5(4):259-76.
5. Glock et al. Experimental techniques of aorto-ostomy anastomosis by gluing (BioGlue - CryoLife). *Laparoscopic Aortic/Thoracic Surgery for Occlusive Disease and Aneurysms*. Jan. 28, 2000.
6. Hasse J, et al. Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. *European J Cardio-Thorac Surg* 2001;19(1):4-9.
7. Kumar A, et al. Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (In Press) *J Clinical Neuroscience*.
8. Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transphenoidal procedures. *J Clinical Neuroscience* 2003;10(1):92-95.
9. Lang H, et al. Use of BioGlue in nephron sparing surgery for renal cell carcinomas (RCC). A pilot study. Presented at the XVIIIth European Association of Urology on 13 March 2003 (Abstract).
10. LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. Presented at the 34th Association of Academic Surgery 2000.
11. LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. Presented at Symposium for Cardiac and General Thoracic Surgery Teams. (Abstract) 2000.
12. LeMaire SA, et al. BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. *Ann Thorac Surg* 2002;73(5):1500-5.
13. LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
14. LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
15. Passage J, et al. BioGlue surgical adhesive-an appropriate tool for indications in cardiac surgery. *Ann Thorac Surg* 2002;74:432-7.
16. Potaris K, et al. Experience with an albumin-glutaraldehyde tissue adhesive in sealing air leaks after bullectomy. *Heart Surgery Forum* 2003;6(5):429-433.
17. Potaris K, et al. Preliminary results with the use of an albumin-glutaraldehyde tissue adhesive in lung surgery. *Medical Science Monitor* 2003;9(7):179-83.
18. Raanani E, et al. Use of "BioGlue" in aortic surgical repair. *Ann Thorac Surg* 2001;72:638-40.
19. Sen A, et al. Cerebrospinal fluid leak rate after the use of BioGlue in transabyrinthine vestibular schwannoma surgery: A prospective study. *Otol Neurotol* 2006;27:102-5.
20. Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. *J Thorac Cardiovasc Surg* 2006;132:105-12.
21. Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. *Ann Thorac Surg* 2002;73:707-13.

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- Bavaria JE, et al. Advances in the treatment of acute type A dissection: an integrated approach. Ann Thorac Surg 2002;74:S1848-52.
- Bavaria JE, et al. Aortic arch dissection: current perspectives. J Penn. Assoc. Thorac Surg 1999;7-8.
- Coselli JS, et al. Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. J Am Coll Surg 2003;197(2):243-253.
- Coselli JS, et al. Thoracic aortic anastomoses. Operative Techniques in Thoracic and Cardiovascular Surg 2000;5(4):259-76.
- Glock et al. Experimental techniques of aortic prosthetic anastomosis by gluing (BioGlue – Cryolife). Laparoscopic Aorticliac Surgery for Occlusive Disease and Aneurysms. Jan. 28, 2000.
- Hasse J, et al. Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. European J Cardio-Thorac Surg 2001;19(1):4-9.
- Kumar A, et al. Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (In Press) J Clinical Neuroscience.
- Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transphenoidal procedures. J Clinical Neuroscience 2003;10(1):92-95.
- Lang H, et al. Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). A pilot study. Presented at the XVIIIth European Association of Urology on 13 March 2003 (Abstract).
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes sciatic nerve injury and diaphragmatic paralysis. Presented at the 34th Association of Academic Surgery 2000.
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. Presented at Symposium for Cardiac and General Thoracic Surgery Teams. (Abstract) 2000.
- LeMaire SA, et al. BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. Ann Thorac Surg 2002;73(5):1500-5.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue surgical adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
- Passage J, et al. BioGlue surgical adhesive-an appraisal of its indications in cardiac surgery. Ann Thorac Surg 2002;74:432-7.
- Potaris K, et al. Experience with an albumin-glutaraldehyde tissue adhesive in sealing air leaks after bullectomy. Heart Surgery Forum 2003;6(5):429-433.
- Potaris K, et al. Preliminary results with the use of an albumin-glutaraldehyde tissue adhesive in lung surgery. Medical Science Monitor 2003;9(7):179-83.
- Raanani E, et al. Use of "BioGlue" in aortic surgical repair. Ann Thorac Surg 2001;72:638-40.
- Sen A, et al. Cerebrospinal fluid leak rate after the use of BioGlue in translabrynthine vestibular schwannoma surgery: A prospective study. Otol Neurotol 2006;27:102-5.
- Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. J Thorac Cardiovasc Surg 2006;132:105-12.
- Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. Ann Thorac Surg 2002;73:707-13.

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- Bavaria JE, et al. Advances in the treatment of acute type A dissection: an integrated approach. Ann Thorac Surg 2002;74:S1848-52.
- Bavaria JE. Aortic aortic dissection: current perspectives. J Penn. Assoc. Thorac Surg 1999;7-8.
- Coselli JS, et al. Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. J Am Coll Surg 2003;197(2):243-253.
- Coselli JS, et al. Thoracic aortic anastomoses. Operative Techniques in Thoracic and Cardiovascular Surg 2000;5(4):259-76.
- Glock et al. Experimental techniques of aorto-prosthetic anastomosis by gluing (BioGlue - CryoLife). Laparoscopic Aortic/iliac Surgery for Occlusive Disease and Aneurysms. Jan. 28, 2000.
- Hasse J, et al. Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. European J Cardio-Thorac Surg 2001;19(1):4-9.
- Kumar A, et al. Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (In Press) J Clinical Neuroscience.
- Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transphenoidal procedures. J Clinical Neuroscience 2003;10(1):92-95.
- Lang H, et al. Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). A pilot study. Presented at the XVIIIth European Association of Urology on 13 March 2003 (Abstract).
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. Presented at the 34th Association of Academic Surgery 2000.
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. Presented at Symposium for Cardiac and General Thoracic Surgery Teams. (Abstract) 2000.
- LeMaire SA, et al. BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. Ann Thorac Surg 2002;73(5):1500-5.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
- Passage J, et al. BioGlue surgical adhesive-an appraisal of its indications in cardiac surgery. Ann Thorac Surg 2002;74:432-7.
- Polaris K, et al. Experience with an albumin-glutaraldehyde tissue adhesive in sealing air leaks after bullectomy. Heart Surgery Forum 2003;6(5):429-433.
- Polaris K, et al. Preliminary results with the use of an albumin-glutaraldehyde tissue adhesive in lung surgery. Medical Science Monitor 2003;9(7):179-83.
- Raanan E, et al. Use of "BioGlue" in aortic surgical repair. Ann Thorac Surg 2001;72:638-40.
- Sen A, et al. Cerebrospinal fluid leak rate after the use of BioGlue in translabrynthine vestibular schwannoma surgery: A prospective study. Otol Neurotol 2006;27:102-5.
- Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. J Thorac Cardiovasc Surg 2006;132:105-12.
- Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. Ann Thorac Surg 2002;73:707-13.

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1. Bavaria JE, et al. Advances in the treatment of acute type A dissection: an integrated approach. *Ann Thorac Surg* 2002;74:S1848-52.
2. Bavaria JE. Acute aortic dissection: current perspectives. *J Penn. Assoc. Thorac Surg* 1999;7-8.
3. Coselli JS, et al. Prospective randomized study of a protein-based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. *J Am Coll Surg* 2003;197(2):243-253.
4. Coselli JS, et al. Thoracic aortic anastomoses. Operative Techniques in Thoracic and Cardiovascular Surg 2000;5(4):259-76.
5. Glock et al. Experimental techniques of aorto-prosthetic anastomosis by gluing (BioGlue - CryoLife). *Laparoscopic Aortic/ail Surgery for Occlusive Disease and Aneurysms*. Jan. 28, 2000.
6. Hasse J, et al. Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. *European J Cardio-Thorac Surg* 2001;19(1):4-9.
7. Kumar A, et al. Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (In Press) *J Clinical Neurosciences*.
8. Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transphenoidal procedures. *J Clinical Neuroscience* 2003;10(11):92-95.
9. Lang H, et al. Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). A pilot study. Presented at the XVIIIth European Association of Urology on 13 March 2003 (Abstract).
10. LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. Presented at the 34th Association of Academic Surgery 2000.
11. LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. Presented at Symposium for Cardiac and General Thoracic Surgery Teams. (Abstract) 2000.
12. LeMaire SA, et al. BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. *Ann Thorac Surg* 2002;73(5):1500-5.
13. LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed tissues during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
14. LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
15. Passage J, et al. BioGlue surgical adhesive: an appraisal of its indications in cardiac surgery. *Ann Thorac Surg* 2002;74:432-7.
16. Polaris K, et al. Experience with an albumin-glutaraldehyde tissue adhesive in sealing air leaks after bullectomy. *Heart Surgery Forum* 2003;6(5):429-433.
17. Polaris K, et al. Preliminary results with the use of an albumin-glutaraldehyde tissue adhesive in lung surgery. *Medical Science Monitor* 2003;9(7):179-83.
18. Raanani E, et al. Use of "BioGlue" in aortic surgical repair. *Ann Thorac Surg* 2001;72:638-40.
19. Sen A, et al. Cerebrosplinal fluid leak rate after the use of BioGlue in translabrynthine vestibular schwannoma surgery: A prospective study. *Otol Neurotol* 2006;27:102-5.
20. Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. *J Thorac Cardiovasc Surg* 2006;132:105-12.
21. Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. *Ann Thorac Surg* 2002;73:707-13.

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- Bavaria JE, et al. Advances in the treatment of acute type A dissection: an integrated approach. *Ann Thorac Surg* 2002;74(S1):84-82.
- Bavaria JE. Acute aortic dissection: current perspectives. *J Penn. Assoc. Thorac Surg* 1999;7:9-8
- Coselli JS, et al. Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. *J Am Coll Surg* 2003;197(2):243-253-
- Coselli JS, et al. Thoracic aortic anastomoses. *Operative Techniques in Thoracic and Cardiovascular Surg* 2000;5(4):259-76.
- Glock et al. Experimental techniques of aorto-prosthetic anastomosis by gluing (BioGlue – CryoLife). *Laparoscopic Aortic Surgery for Occlusive Disease and Aneurysms*. 28 Jan., 2000.
- Hasse J, et al. Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. *European J Cardio-Thorac Surg* 2001;19(1):4-9.
- Kumar A, et al. Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (In Press) *J Clinical Neuroscience*
- Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transphenoidal procedures. *J Clinical Neuroscience* 2003;10(1):92-95.
- Lang H, et al. Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). *En pilotstudie*. *Presentatades vid 18th European Association of Urology* den 13 mars, 2003. (utdrag).
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. *Presentatades vid 34th Association of Academic Surgery* 2000.
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. *Presentatades vid Symposium for Cardiac and General Thoracic Surgery Teams*. (utdrag) 2000.
- LeMaire SA, et al. BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. *Ann Thorac Surg* 2002;73(5):1500-5.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. *Presentatades vid 35th meeting for Association of Academic Surgery*. (utdrag) 2001.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. *Presentatades vid 35th meeting for Association of Academic Surgery*. (utdrag) 2001.
- Passage J, et al. BioGlue surgical adhesive-an appraisal of its indications in cardiac surgery. *Ann Thorac Surg* 2002;74:432-7.
- Polaris K, et al. Experience with an albumin-glutaraldehyde tissue adhesive in sealing air leaks after bullectomy. *Heart Surgery Forum* 2003;6(5):429-433.
- Polaris K, et al. Preliminary results with the use of an albumin-glutaraldehyde tissue adhesive in lung surgery. *Medical Science Monitor* 2003;9(7):179-83.
- Raanani E, et al. Use of "BioGlue" in aortic surgical repair. *Ann Thorac Surg* 2001;72:638-40.
- Sen A, et al. Cerebrosplinal fluid leak rate after the use of BioGlue in translabrynthine vestibular schwannoma surgery: A prospective study. *Otol Neurotol* 2006;27:102-5.
- Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. *J Thorac Cardiovasc Surg* 2006;132:105-12.
- Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. *Ann Thorac Surg* 2002;73:707-13.

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- Bavaria JE, et al. Advances in the treatment of acute type a aortic dissection: an integrated approach. Ann Thorac Surg 2002;74:1848-52.
- Bavaria JE. Acute aortic dissection: current perspectives. J Penn. Assoc. Thorac Surg 1999;7-8.
- Coselli JS, et al. Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. J Am Coll Surg 2003;197(2):243-253.
- Coselli JS, et al. Thoracic aortic anastomoses. Operative Techniques in Thoracic and Cardiovascular Surg 2003;5(4):259-76.
- Glock et al. Experimental techniques of aortic-posterior anastomosis by gluing (BioGlue – CryoLife). Laparoscopic Aortic/Urg for Occlusive Disease and Aneurysms. Jan. 28, 2000.
- Hasse J, et al. Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. European J Cardio-Thorac Surg 2001;19(1):4-9.
- Kumar A, et al. Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (In Press) J Clinical Neuroscience.
- Kumar A, et al. Reconstruction of the sellar floor using BioGlue in transphenoidal procedures. J Clinical Neuroscience 2003;10(1):92-95.
- Lang H, et al. Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). A pilot study. Presented at the XVIIIth European Association of Urology on 13 March 2003 (Abstract).
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. Presented at the 34th Association of Academic Surgery 2000.
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged diaphragmatic paralysis due to sinoatrial node degeneration. Presented at Symposium for Cardiac and General Thoracic Surgery Teams. (Abstract) 2000.
- LeMaire SA, et al. BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. Ann Thorac Surg 2002;73(5):1500-5.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
- Passage J, et al. BioGlue surgical adhesive-an appraisal of its indications in cardiac surgery. Ann Thorac Surg 2002;74:432-7.
- Polaris K, et al. Experience with an albumin-glutaraldehyde tissue adhesive in sealing air leaks after bullectomy. Heart Surgery Forum 2003;6(5):429-433.
- Polaris K, et al. Preliminary results with the use of an albumin-glutaraldehyde tissue adhesive in lung surgery. Medical Science Monitor 2003;9(7):179-83.
- Raanan E, et al. Use of "BioGlue" in aortic surgical repair. Ann Thorac Surg 2001;72:638-40.
- Sen A, et al. Cerebrosplinal fluid leak rate after the use of BioGlue in translabrynthine vestibular schwannoma surgery: A prospective study. Otol Neurotol 2006;27:102-5.
- Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. J Thorac Cardiovasc Surg 2006;132:105-12.
- Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. Ann Thorac Surg 2002;73:707-13.

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- Bavaria JE, et al. Advances in the treatment of acute type A dissection: an integrated approach. Ann Thorac Surg 2002;74:S1846-52.
- Bavaria JE. Acute aortic dissection: current perspectives. J Penn. Assoc. Thorac Surg 1999;7-8.
- Coselli JS, et al. Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. J Am Coll Surg 2003;197(2):243-253.
- Coselli JS, et al. Thoracic aortic anastomoses. Operative Techniques in Thoracic and Cardiovascular Surgery 2000;5(4):259-76.
- Glock et al. Experimental techniques of aorto-prosthetic anastomosis by gluing (BioGlue - CryoLife). Laparoscopic Aortic/iaic Surgery for Occlusive Disease and Aneurysms. Ocak 28, 2000.
- Hasse J, et al. Experimental use of an albumin-glyceraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. European J Cardio-Thorac Surg 2001;19(1):4-9.
- Kumar A, et al. Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (In Press) J Clinical Neuroscience.
- Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transphenoidal procedures. J Clinical Neuroscience 2003;10(1):92-95.
- Lang H, et al. Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). A pilot study, 13 Mart 2003 tarihinde XVIII. Avrupa Üroloji Birliğinde sunulmuştur (Özet).
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. 34. Akademik Cerrahi Birliği 2000 kapaşında sunulmuştur.
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. Kalp ve Genel Göğüs Cerrahisi Ekipleri Sempozyumunda sunulmuştur. (Özet) 2000.
- LeMaire SA, et al. BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. Ann Thorac Surg 2002;73(5):1500-5.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. 35. Akademik Cerrahi Birliği Toplantısında sunulmuştur. (Özet) 2001.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. 35. Akademik Cerrahi Birliği Toplantısında sunulmuştur. (Özet) 2001.
- Passage J, et al. BioGlue surgical adhesive-an appraisal of its indications in cardiac surgery. Ann Thorac Surg 2002;74:432-7.
- Polaris K, et al. Experience with an albumin-glyceraldehyde tissue adhesive in sealing air leaks after bullectomy. Heart Surg Forum 2003;6(5):429-433.
- Polaris K, et al. Preliminary results with the use of an albumin-glyceraldehyde tissue adhesive in lung surgery. Medical Science Monitor 2003;9(7):179-83.
- Raamani E, et al. Use of "BioGlue" in aortic surgical repair. Ann Thorac Surg 2001;72:638-40.
- Sen A, et al. Cerebrospinal fluid leak rate after the use of BioGlue in translabrynthine vestibular schwannoma surgery: A prospective study. Otol Neurotol 2006;27:102-5.
- Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. J Thorac Cardiovasc Surg 2006;132:105-12.
- Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. Ann Thorac Surg 2002;73:707-13.

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1. Bavaria JE, et al. Advances in the treatment of acute type A dissection: an integrated approach. *Ann Thorac Surg* 2002;74:S1848-52.
2. Bavaria JE. Acute aortic dissection: current perspectives. *J Penn. Assoc. Thorac Surg* 1999;7:8.
3. Coselli JS, et al. Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. *J Am Coll Surg* 2003;197(2):243-253.
4. Coselli JS, et al. Thoracic aortic anastomoses. *Operative Techniques in Thoracic and Cardiovascular Surg* 2000;5(4):259-76.
5. Glock et al. Experimental techniques of aorto-prosthetic anastomosis by gluing (BioGlue - CryoLife). *Laparoscopic Aortic/iliac Surgery for Occlusive Disease and Aneurysms*. Jan. 28, 2000.
6. Hasse J, et al. Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. *European J Cardio-Thorac Surg* 2001;19(1):4-9.
7. Kumar A, et al. Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (In Press) *J Clinical Neuroscience*.
8. Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transphenoidal procedures. *J Clinical Neuroscience* 2003;10(1):92-95.
9. Lang H, et al. Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). A pilot study. Presented at the XVIIIth European Association of Urology on 13 March 2003 (Abstract).
10. LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. Presented at the 34th Association of Academic Surgery 2000.
11. LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. Presented at Symposium for Cardiac and General Thoracic Surgery Teams. (Abstract) 2000.
12. LeMaire SA, et al. BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. *Ann Thorac Surg* 2002;73(5):1500-5.
13. LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
14. LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
15. Passage J, et al. BioGlue surgical adhesive-an appraisal of its indications in cardiac surgery. *Ann Thorac Surg* 2002;74:432-7.
16. Polaris K, et al. Experience with an albumin-glutaraldehyde tissue adhesive in sealing air leaks after bullectomy. *Heart Surgery Forum* 2003;6(5):429-433.
17. Polaris K, et al. Preliminary results with the use of an albumin-glutaraldehyde tissue adhesive in lung surgery. *Medical Science Monitor* 2003;9(7):179-83.
18. Raanani E, et al. Use of "BioGlue" in aortic surgical repair. *Ann Thorac Surg* 2001;72:638-40.
19. Sen A, et al. Cerebrospinal fluid leak rate after the use of BioGlue in translabirithine vestibular schwannoma surgery: A prospective study. *Otol Neurotol* 2006;27:102-5.
20. Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. *J Thorac Cardiovasc Surg* 2006;132:105-12.
21. Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. *Ann Thorac Surg* 2002;73:707-13.

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1. Bavaria JE, et al. Advances in the treatment of acute type A dissection: an integrated approach. Ann Thorac Surg 2002;74:S184-82.
2. Bavaria JE. Acute aortic dissection: current perspectives. J Penn. Assoc. Thorac Surg 1999;7-8.
3. Coselli JS, et al. Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. J Am Coll Surg 2003;197(2):243-253.
4. Coselli JS, et al. Thoracic aortic anastomoses. Operative Techniques in Thoracic and Cardiovascular Surg 2000;5(4):259-76.
5. Glock et al. Experimental techniques of aorto-prosthetic anastomosis by gluing (BioGlue - CryoLife). Laparoscopic Aortic/Inf Surgery for Occlusive Disease and Aneurysms. Ozak 28, 2000.
6. Hassa J, et al. Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. European J Cardio-Thorac Surg 2001;19(1):4-9.
7. Kumar A, et al. Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (In Press) J Clinical Neuroscience.
8. Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transphenoidal procedures. J Clinical Neuroscience 2003;10(1):92-95.
9. Lang H, et al. Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). A pilot study. 13 Mar 2003 Larihinde XVIII. Avrupa Uroloji Birliđinde sunulmuştur (Özet).
10. LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phase nerve injury and diaphragmatic paralysis. 34. Akademik Cerrahi Birliđi 2000 kaptamında sunulmuştur.
11. LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. Kalp ve Genel Göğüs Cerrahisi Ekipleri Sempozyumu'nda sunulmuştur. (Özet) 2000.
12. LeMaire SA, et al. BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. Ann Thorac Surg 2002;73(5):1500-5.
13. LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. 35. Akademik Cerrahi Birliđi Toplantısında sunulmuştur. (Özet) 2001.
14. LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. 35. Akademik Cerrahi Birliđi Toplantısında sunulmuştur. (Özet) 2001.
15. Passage J, et al. BioGlue surgical adhesive-an appraisal of its indications in cardiac surgery. Ann Thorac Surg 2002;74:432-7.
16. Polaris K, et al. Experience with an albumin-glutaraldehyde tissue adhesive in sealing air leaks after bullectomy. Heart Surgery Forum 2003;6(5):429-433.
17. Polaris K, et al. Preliminary results with the use of an albumin-glutaraldehyde tissue adhesive in lung surgery. Medical Science Monitor 2003;9(7):179-83.
18. Raanan E, et al. Use of "BioGlue" in aortic surgical repair. Ann Thorac Surg 2001;72:638-40.
19. Sen A, et al. Cerebrospinal fluid leak rate after the use of BioGlue in transabyrinthine vestibular schwannoma surgery: A prospective study. Otol Neurotol 2006;27:102-5.
20. Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. J Thorac Cardiovasc Surg 2006;132:105-12.
21. Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. Ann Thorac Surg 2002;73:707-13.

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1. Bavaria JE, et al. Advances in the treatment of acute type A dissection: an integrated approach. Ann Thorac Surg 2002;74(1):848-52.
2. Bavaria JE. Acute aortic dissection: current perspectives. J Penn. Assoc. Thorac Surg 1999;7-8.
3. Coselli JS, et al. Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. J Am Coll Surg 2003;197(2):243-253.
4. Coselli JS, et al. Thoracic aortic anastomoses. Operative Techniques in Thoracic and Cardiovascular Surg 2000;5(4):259-76.
5. Glock et al. Experimental techniques of aorto-prosthetic anastomosis by gluing (BioGlue - CryoLife). Laparoscopic Aortic/Infra Surgery for Occlusive Disease and Aneurysms. Ocak 28, 2000.
6. Hesse J, et al. Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. European J Cardio-Thorac Surg 2001;19(1):4-9.
7. Kumar A, et al. Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (In Press) J Clinical Neuroscience.
8. Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transsphenoidal procedures. J Clinical Neuroscience 2003;10(1):92-95.
9. Lang H, et al. Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). A pilot study. 13 Mart 2003 tarinhde XVIII. Avrupa Uroloji Birliđinde sunulmuřtur (Özet).
10. LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. 34. Akademik Cerrahi Birliđi 2000 kapsamında sunulmuřtur.
11. LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. Kalp ve Genel Gögüs Cerrahisi Ekipleri Sempozyumu'nda sunulmuřtur. (Özet) 2000.
12. LeMaire SA, et al. BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. Ann Thorac Surg 2002;73(5):1500-5.
13. LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. 35. Akademik Cerrahi Birliđi Toplantısında sunulmuřtur. (Özet) 2001.
14. LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. 35. Akademik Cerrahi Birliđi Toplantısında sunulmuřtur. (Özet) 2001.
15. Passage J, et al. BioGlue surgical adhesive-an appraisal of its indications in cardiac surgery. Ann Thorac Surg 2002;74:432-7.
16. Potaris K, et al. Experience with an albumin-glutaraldehyde tissue adhesive in sealing air leaks after bullectomy. Heart Surgery Forum 2003;6(5):429-433.
17. Potaris K, et al. Preliminary results with the use of an albumin-glutaraldehyde tissue adhesive in lung surgery. Medical Science Monitor 2003;9(7):179-83.
18. Raanani E, et al. Use of "BioGlue" in aortic surgical repair. Ann Thorac Surg 2001;72:638-40.
19. Sen A, et al. Cerebrospinal fluid leak rate after the use of BioGlue in translabrynthine vestibular schwannoma surgery: A prospective study. Otol Neurotol 2006;27:102-5.
20. Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. J Thorac Cardiovasc Surg 2006;132:105-12.
21. Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. Ann Thorac Surg 2002;73:707-13.

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1. Bavaria JE, et al. Advances in the treatment of acute type A dissection: an integrated approach. *Ann Thorac Surg* 2002;74:1848-52.
2. Bavaria JE. Acute aortic dissection: current perspectives. *J Penn. Assoc. Thorac Surg* 1999;7-8.
3. Coselli JS, et al. Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. *J Am Coll Surg* 2003;197(2):243-253.
4. Coselli JS, et al. Thoracic aortic anastomoses. Operative Techniques in Thoracic and Cardiovascular Surg 2000;5(4):259-76.
5. Glock et al. Experimental techniques of aorto-prosthetic anastomosis by gluing (BioGlue - CryoLife). *Laparoscopic Aortic/Clavicular Surgery for Occlusive Disease and Aneurysms*. *Oak* 28, 2000.
6. Hassa J, et al. Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. *European J Cardio-Thorac Surg* 2001;19(1):4-9.
7. Kumar A, et al. Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (In Press) *J Clinical Neuroscience*.
8. Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transphenoidal procedures. *J Clinical Neuroscience* 2003;10(1):92-95.
9. Lang H, et al. Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). A pilot study. 13 Mar 2003 Iarhinde XVIII. *Avrupa Urologi Birliğinde sunulmuşdur* (Özet).
10. LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. 34. *Akademik Cerrahi Birliği 2000 kampsında sunulmuşdur*.
11. LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. *Kalp ve Genel Göğüs Cerrahisi Etkileri Sempozyumu'nda sunulmuşdur*. (Özet) 2000.
12. LeMaire SA, et al. BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. *Ann Thorac Surg* 2002;73(5):1500-5.
13. LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. 35. *Akademik Cerrahi Birliği Toplantısında sunulmuşdur*. (Özet) 2001.
14. LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. 35. *Akademik Cerrahi Birliği Toplantısında sunulmuşdur*. (Özet) 2001.
15. Passage J, et al. BioGlue surgical adhesive-an appraisal of its indications in cardiac surgery. *Ann Thorac Surg* 2002;74:432-7.
16. Polaris K, et al. Experience with an albumin-glutaraldehyde tissue adhesive in sealing air leaks after bullectomy. *Heart Surgery Forum* 2003;6(5):429-433.
17. Polaris K, et al. Preliminary results with the use of an albumin-glutaraldehyde tissue adhesive in lung surgery. *Medical Science Monitor* 2003;9(7):179-83.
18. Raanan E, et al. Use of "BioGlue" in aortic surgical repair. *Ann Thorac Surg* 2001;72:638-40.
19. Sen A, et al. Cerebrospinal fluid leak rate after the use of BioGlue in translabrynthine vestibular schwannoma surgery: A prospective study. *Otol Neurotol* 2006;27:102-5.
20. Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. *J Thorac Cardiovasc Surg* 2006;132:105-12.
21. Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. *Ann Thorac Surg* 2002;73:707-13.

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- Bavaria JE, et al. Advances in the treatment of acute type A dissection: an integrated approach. Ann Thorac Surg 2002;74:S1848-52.
- Bavaria JE. Acute aortic dissection: current perspectives. J Penn. Assoc. Thorac Surg 1999;7-8.
- Coselli JS, et al. Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. J Am Coll Surg 2003;197(2):243-253.
- Coselli JS, et al. Thoracic aortic anastomoses. Operative Techniques in Thoracic and Cardiovascular Surg 2000;5(4):258-71.
- Glock et al. Experimental techniques of aortic prosthetic anastomosis by gluing (BioGlue – CryoLife). Laparoscopic Aorticliac Surgery for Occlusive Disease and Aneurysms. Ooak 28, 2000.
- Hasse J, et al. Experimental use of an albumin-glyceraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. European J Cardio-Thorac Surg 2001;19(1):4-9.
- Kumar A, et al. Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (In Press) J Clinical Neuroscience.
- Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transsphenoidal procedures. J Clinical Neuroscience 2005;10(1):92-95.
- Lang H, et al. Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). A pilot study. 13 Mart 2003 tarhinde XVIII. Avrupa Uroloji Birliği'nde sunulmuştur (Özet).
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. 34. Akademik Cerrahi Birliği 2000 kampsında sunulmuştur.
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. Kalp ve Genel Göğüs Cerrahisi Ekipleri Sempozyumu'nda sunulmuştur. (Özet) 2000.
- LeMaire SA, et al. BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. Ann Thorac Surg 2002;74(5):1500-5.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. 35. Akademik Cerrahi Birliği Toplantısında sunulmuştur. (Özet) 2001.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. 35. Akademik Cerrahi Birliği Toplantısında sunulmuştur. (Özet) 2001.
- Passage J, et al. BioGlue surgical adhesive-an appraisal of its indications in cardiac surgery. Ann Thorac Surg 2002;74:432-7.
- Polaris K, et al. Experience with an albumin-glyceraldehyde tissue adhesive in sealing air leaks after bullectomy. Heart Surgery Forum 2003;6(5):423-429.
- Polaris K, et al. Preliminary results with the use of an albumin-glyceraldehyde tissue adhesive in lung surgery. Medical Science Monitor 2003;9(7):179-83.
- Raanani E, et al. Use of "BioGlue" in aortic surgical repair. Ann Thorac Surg 2001;72:638-40.
- Sen A, et al. Cerebrospinal fluid leak rate after the use of BioGlue in transabyrinthine vestibular schwannoma surgery: A prospective study. Otol Neurotol 2006;27:102-5.
- Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating mediastinal air leaks. J Thorac Cardiovasc Surg 2006;132:105-12.
- Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. Ann Thorac Surg 2002;73:707-13.

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- Bavaria JE, et al. Advances in the treatment of acute type A dissection: an integrated approach. *Ann Thorac Surg* 2002;74:51848-52.
- Bavaria JE. Acute aortic dissection: current perspectives. *J Penn Assoc. Thorac Surg* 1999;7-8.
- Coselli JS, et al. Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. *J Am Coll Surg* 2003;197(2):243-253.
- Coselli JS, et al. Thoracic aortic anastomoses. *Operative Techniques in Thoracic and Cardiovascular Surg* 2000;5(4):259-76.
- Glock et al. Experimental techniques of aorto-prosthetic anastomosis by gluing (BioGlue – CryoLife). *Laparoscopic Aortic Surgery for Occlusive Disease and Aneurysms*. Jan 28, 2000.
- Hasse J, et al. Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. *European J Cardio-Thorac Surg* 2001;19(1):4-9.
- Kumar A, et al. Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (In Press) *J Clinical Neuroscience*.
- Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transphenoidal procedures. *J Clinical Neuroscience* 2003;10(1):92-95.
- Lang H, et al. Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). A pilot study. Presented at the XXVIII European Association of Urology on 13 March 2003 (Abstract).
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. Presented at the 34th Association of Academic Surgery 2000.
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. Presented at Symposium for Cardiac and General Thoracic Surgery Teams. (Abstract) 2000.
- LeMaire SA, et al. BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. *Ann Thorac Surg* 2002;73(5):1500-5.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
- Passage J, et al. BioGlue surgical adhesive—an appraisal of its indications in cardiac surgery. *Ann Thorac Surg* 2002;74:432-7.
- Potaris K, et al. Experience with an albumin-glutaraldehyde tissue adhesive in sealing air leaks after bullectomy. *Heart Surgery Forum* 2003;6(5):429-433.
- Potaris K, et al. Preliminary results with the use of an albumin-glutaraldehyde tissue adhesive in lung surgery. *Chest* 2003;123(1):179-83.
- Raanani E, et al. Use of "BioGlue" in aortic surgical repair. *Ann Thorac Surg* 2001;72:638-40.
- Sen A, et al. Cerebrospinal fluid leak rate after the use of BioGlue in transabyrinthine vestibular schwannoma surgery: A prospective study. *Otol Neurotol* 2006;27:102-5.
- Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. *J Thorac Cardiovasc Surg* 2006;132:105-12.
- Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. *Ann Thorac Surg* 2002;73:707-13.

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Representative for vigilance reports (except for European Union Countries):

CryoLife, Inc. • 1655 Roberts Blvd., NW • Kennesaw, Georgia 30144, USA
Phone: +1 (770) 419-3355 • FAX +1 (770) 590-3753 • E-mail: usbioglu@cryolife.com

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CryoLife, Inc. • 1655 Roberts Boulevard, NW • Kennesaw, Georgia 30144, USA