### Delivery Tip Extension

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= Sterilized Using Ethylene Oxide	= Authorised Representative in the European Community
= Note made with natural rubber latex	= Manufacturer
= Non-Pyrogenic	= Use-by date
= Do Not Re-use	= Catalogue Number
= Do Not Use if Package is Damaged	= Batch Code
= Consult Instructions for Use	= Prescription Device
= Caution	

In the United States of America Federal law restricts a prescription device to sale by or on the order of a physician.  
In all other regions of the world, this device should not be used except on the order of a physician.

Read the *Instructions for Use* prior to using this product.

## **DEVICE DESCRIPTION**

The delivery tip extension (DTE) is intended to be used with BioGlue® Surgical Adhesive (BioGlue) syringe standard applicator tips. The DTEs, available in a variety of lengths to accommodate different surgical needs, slip over the standard applicator tips to deliver BioGlue for use as stated in the BioGlue product *Instructions for Use* (IFU). The DTE is made from biocompatible, USP Class VI compliant polymer resins.

## **INTENDED USE**

The DTE is intended for use with BioGlue syringe standard applicator tips to deliver BioGlue to the surgical site.

Please refer to the BioGlue IFU for additional information such as indications for use, contraindications, warnings, and precautions.

## **PACKAGING AND STORAGE**

The DTEs are double pouched (inner and outer pouch: Tyvek®/Poly) with 2 DTEs per pouch and 10 pouches per box and are supplied sterile for single-patient use only. Discard damaged or unused product components from open or compromised packaging.

Product should be stored at room temperature.

## **DIRECTIONS FOR USE**

### BioGlue Preparation

Please refer to and follow the detailed instructions in the BioGlue Syringe IFU for **Device Preparation** Steps 1-5, **Site Preparation**, and **Syringe Air Space Removal** Steps 1 and 2.

**NOTE:** Continue to hold the syringe upright during the entire assembly of the BioGlue delivery system and syringe air space removal.

### **Attaching the Delivery Tip Extension**

1. Within the sterile field, carefully remove a DTE from the packaging.

**NOTE:** The DTE is sold separately.

2. Insert the flared end of the DTE over the applicator tip and push until fully seated. The DTE is fully seated when the tapered end of the standard applicator tip contacts the inner tapered portion of the DTE, as shown in Figure 1.


**Figure 1 – Fully seated**

**Figure 2 – Not fully seated**

**CAUTION:** Do not compress the plunger while attaching the DTE to the applicator tip. Compressing the plunger can be avoided by: (1) Removing the plunger while attaching the DTE and reinserting the plunger after the DTE is attached,

or (2) Taking care not to manually compress the plunger or compress it against any surface.

**Figure 3 – Wrong**

**Figure 4 – Right**

### **BioGlue Applicator Tip and Delivery Tip Extension Preparation and Priming**

1. Each applicator tip/DTE assembly must be primed prior to BioGlue application. Priming ensures that the BioGlue solutions are properly mixed. The surgeon should compress the plunger and expel a narrow ribbon of BioGlue approximately 3 cm long onto a sterile disposable surface (e.g., sponge, gauze or towel).
2. The surgeon should examine the material expelled during priming and ensure that it is of uniform light yellow to amber color and that it is free from air bubbles. If this material looks colorless or contains bubbles, repeat the prime as outlined in Step 1 until the device delivers a uniform liquid with no bubbles.

**CAUTION:** Avoid direct contact with material expelled during priming.

3. When the applicator tip/DTE assembly has been properly primed, proceed immediately to application.

**CAUTION:** BioGlue polymerizes very quickly. The surgeon must apply BioGlue immediately after priming. Pausing between priming and application can cause polymerization of BioGlue within the applicator tip/DTE assembly. Should this occur, replace the obstructed applicator tip/DTE assembly with a new applicator tip and attach another DTE; then repeat the steps for priming. Do not continue to apply pressure to the plunger once the tip has occluded.

**NOTE:** The DTE is intended to be used to facilitate delivery of BioGlue and is not intended to be used as a surgical tool (i.e. pressing, lifting, and/or holding tissues).

### **PRODUCT INFORMATION DISCLOSURE**

Handling and storage of this device by the user, as well as the factors related to the patient, the patient's diagnosis, treatment, surgical procedures, and other matters beyond manufacturer's control, may directly or indirectly affect this device and the results obtained from its use. This device should not be used except on the order of a physician.

### **DISCLAIMER OF WARRANTIES; LIMITS OF LIABILITY**

CRYOLIFE DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES WITH RESPECT TO THIS DELIVERY TIP EXTENSION, INCLUDING BUT NOT LIMITED TO THE EXPRESS AND IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR

PURPOSE. IN NO EVENT SHALL CRYOLIFE BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES. In the event that such disclaimer is found invalid or unenforceable for any reason: (i) any action for breach of warranty must be commenced within one year after any such claim or cause of action accrued and (ii) the remedy for any such breach is limited to the replacement of the product. Prices, specifications, and availability are subject to change without notice.

<b>CE0088</b> = CE Marking of Conformity	<b>EC REP</b> = Authorised Representative in the European Community
<b>STERILE EO</b> = Sterilized Using Ethylene Oxide	 = Manufacturer
<b>LATEX</b> = Not made with natural rubber latex	 = Use-by date
 = Non-Pyrogenic	<b>REF</b> = Catalogue Number
 = Do Not Re-use	<b>LOT</b> = Batch Code
 = Do Not Use if Package is Damaged	<b>RxONLY</b> = Prescription Device
 = Consult Instructions for Use	
 = Caution	

*In the United States of America Federal law restricts a prescription device to sale by or on the order of a physician.*

*In all other regions of the world, this device should not be used except on the order of a physician.*

**Representative for vigilance reports (except for European Union countries):**

 CryoLife, Inc.  
 1655 Roberts Boulevard, NW  
 Kennesaw, Georgia 30144 • USA  
 1-800-438-8285 (United States and Canada)  
 +1-770-419-3355 (International)  
 +1-770-590-3753 (Fax)  
[www.cryolife.com](http://www.cryolife.com) (Website)  
[usbioglue@cryolife.com](mailto:usbioglue@cryolife.com) (E-mail)

**EC REP** CryoLife Europa, Ltd. • Bramley House • The Guildway • Old Portsmouth Road • Guildford • Surrey GU3 1LR • United Kingdom  
 Phone: +44 (0) 1483 441030  
 Fax: +44 (0) 1483 452860  
 E-mail: [Europa@cryolife.com](mailto:Europa@cryolife.com)

Manufactured by: CryoLife, Inc.

© 2015 CryoLife, Inc. All rights reserved. The snowflake design, CryoLife, and BioGlue are trademarks of CryoLife, Inc. Tyvek is a registered trademark of E.I. du Pont de Nemours and Company.

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Üretici: CryoLife, Inc.