PhotoFix®

Decellularized Bovine Pericardium

QTY	1
-----	---

Quantity

INSTRUCTIONS FOR USE

Caution: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

PhotoFix® Decellularized Bovine Pericardium (PhotoFix) is prepared from bovine pericardium, which is stabilized using a dye-mediated photooxidation process, processed using ethylene oxide and sterilized using aseptic processing techniques. No aldehyde chemistry is used during any phase of pericardium manufacturing including tissue fixation or sterilization processes.

INDICATIONS

PhotoFix is indicated for the following uses: intracardiac repair, great vessel repair, suture line buttressing, pericardial closure, and vascular repair and reconstruction (for example: the carotid, iliac, femoral, and tibial blood vessels and arteriovenous access revisions).

CONTRAINDICATIONS

PhotoFix is contraindicated for:

- Use in patients that exhibit sensitivity to materials of bovine origin;
- Reconstruction of hernias and valve leaflet repair.

WARNINGS AND PRECAUTIONS

• FOR SINGLE USE ONLY

- PhotoFix is supplied sterile in a sealed container. Do not re-sterilize.
- Do not use the pericardium if either temperature indicator is activated (see **TEMPERATURE INDICATORS**)

section). Immediately contact your Sales Representative to arrange for return or replacement of the tissue.

- **DO NOT USE THE PERICARDIUM** if the device container is damaged, the expiration date has lapsed, or the tamper-evident seal is broken or missing.
- PhotoFix should not be allowed to contact any solutions that contain aldehydes.
- Rinsing the pericardium with antibiotic solution is not recommended. The interactions of tissue with antibiotics have not been tested.

ADVERSE EVENTS

Potential adverse events associated with the use of pericardium include: Abscess, adhesion, calcification/mineralization, complete heart block, death, epicardial inflammatory reactions, infection and rejection.

Adverse events potentially related to the product must be reported promptly to CryoLife Field Assurance at 1-800-438-8285 or 770-419-3355 or by email to fieldassurance@cryolife.com.

INDIVIDUALIZATION OF TREATMENT

Prophylactic antibiotic therapy is necessary for patients that receive an implantable device when undergoing any invasive procedure (including and especially dental procedures).

PATIENT COUNSELING

Patients undergoing any dental procedures should be considered for prophylactic antibiotic therapy.

HOW SUPPLIED

PhotoFix is supplied sterile and non-pyrogenic in a sealed container with 22% buffered ethanol solution. After implantation, the 22% buffered ethanol solution can be discarded according to hospital procedures for non-hazardous materials.

TEMPERATURE INDICATORS

Each shipping carton contains low and high temperature indicators that signify if the tissue was exposed to unacceptable shipping temperatures. If the low temperature indicator is **RED**, the product was exposed to temperatures below the acceptable shipping temperature (**Figure 1**). If the high temperature indicator is **RED**, the product was exposed to temperatures above the acceptable shipping temperature (**Figure 2**).

Figure 1

Figure 2

If either high or low temperature indicators have been activated, **DO NOT USE THE PERICARDIUM**. Immediately contact your Sales Representative to make arrangements for return or replacement.

STORAGE CONDITIONS

PhotoFix should be stored between the temperature limit minimum of 5°C and maximum of 25°C (41°F and 77°F). Refrigeration is not required. Do not place PhotoFix in areas where significant temperature fluctuations may occur.

DIRECTIONS FOR USE

1. The exterior of the container is non-sterile and should not be introduced into the sterile field.
2. Remove the tamper-evident seal by pulling tab adjacent to the seal perforation.
3. Grip the bottle lid and open the container by turning the lid counterclockwise.
4. Remove PhotoFix from the container by grasping the edge of the pericardium with atraumatic forceps. **ATRAUMATIC FORCEPS MUST BE USED WHEN HANDLING PHOTOFIX.** Alternatively, excluding 6x8cm and 8x14cm sizes, the tissue and bottle contents can be poured into a sterile basin.
5. PhotoFix does not require rinsing before implantation.
6. All subsequent handling must be with sterile gloved hands. Thoroughly wash surgical gloves to remove glove powder before touching the pericardium.
7. PhotoFix may be tailored during surgery as needed. Visually examine both sides of PhotoFix. If one side appears smoother, implant the smoother surface so that it faces the blood flow. The pericardium may be sutured, clipped, or stapled into place.
8. Flush the pericardium with sterile physiologic saline frequently during implantation to prevent tissue dehydration.

WARRANTIES

THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO UNITED STATES CUSTOMERS ONLY: ALTHOUGH THE CRYOLIFE, INC. PHOTOFIX® DECELLULARIZED BOVINE PERICARDIUM, HEREAFTER REFERRED TO AS "PRODUCT", HAS BEEN DESIGNED, MANUFACTURED AND TESTED UNDER CAREFULLY CONTROLLED CONDITIONS, DUE TO CONDITIONS AND FACTORS THAT ARE OUTSIDE OF CRYOLIFE, INC.'S SUPERVISION AND CONTROL AFTER SALE OF THE PRODUCT, CRYOLIFE, INC. HEREBY DISCLAIMS ALL WARRANTIES EXPRESSED OR IMPLIED WITH RESPECT TO THE PRODUCT, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. CRYOLIFE, INC. SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY DIRECT, INDIRECT, INCIDENTAL, CONSEQUENTIAL, SPECIAL OR OTHER TYPE OF DAMAGES INCLUDING, WITHOUT LIMITATION, ANY MEDICAL EXPENSES CAUSED BY OR RELATED TO ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER BASED UPON A WARRANTY, CONTRACT, TORT OR OTHER LEGAL THEORY. CRYOLIFE, INC. NEITHER ASSUMES, NOR AUTHORIZES ANY OTHER PARTY TO ASSUME FOR IT, ANY WARRANTY, REPRESENTATION OR OTHER LIABILITY IN CONNECTION WITH THE PRODUCT OR ITS USE.

If any term or portion of this Disclaimer of Warranty and Liability provision is found to be unenforceable, illegal or contrary to applicable law by a court of competent jurisdiction, the remaining portions of this provision shall not be effected and shall remain valid, and any such unenforceability or illegality in jurisdiction shall not invalidate or render unenforceable this Disclaimer of Warranty and Liability provision in any other jurisdiction.

SYMBOL GLOSSARY

RxONLY	=	Prescription only		=	Consult instructions for use
	=	Do not re-sterilize		=	Caution
	=	Sterilized using aseptic processing techniques		=	Temperature limit
	=	Do not re-use		=	Manufacturer
	=	Do not use if package is damaged		=	Catalogue number
	=	Use-by date		=	Batch code

Manufacturer

CryoLife, Inc.
 1655 Roberts Blvd., NW
 Kennesaw, GA 30144
 Telephone: 800-438-8285
 Fax: 770-590-3753

Printed in USA © 2020 CryoLife, Inc., All Rights Reserved.
 PhotoFix is a registered trademark of CryoLife, Inc.

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PhotoFix es una marca comercial registrada de CryoLife, Inc.
© 2020 CryoLife, Inc. Todos los derechos reservados. Impreso en los Estados Unidos.