
ASCYRUS
MEDICAL

AMDS™ Hybrid Prosthesis

Instructions For Use (EN)

Ascyrus Medical GmbH

Grosse Gallusstrasse 16-18
60312 Frankfurt, Germany

0297

A printed IFU will be provided within seven days upon request through any of the contact means listed below:

Phone +1 888.427.9654 • Fax: 770.590.3753 • Email: customerservice.us@artivion.com

1.0 Description

The AMDS Hybrid Prosthesis (AMDS) is constructed of an uncovered Nitinol wire braided stent attached proximally to a PTFE felt graft component made of a PTFE felt tube as shown in Figure 1. The inner diameter of the PTFE felt tube measures 28 mm for AMDS40c and AMDS4030c. The inner diameter of the PTFE felt tube measures 28 mm for AMDS55c and AMDS5540c. The inner diameter of the PTFE felt tube measures 24 mm for AMDS40 and AMDS4030. The inner diameter of the PTFE felt tube measures 32 mm for AMDS55 and AMDS5540. The AMDS is intended for the repair and re-expansion of the intimal flap within the ascending aorta, aortic arch, and into the descending aorta for patients with acute DeBakey Type I aortic dissections and/or intramural hematomas (IMH) undergoing open surgical repair within 0-14 days of diagnosis. Following implantation, in addition to the closure of the proximal part of the false lumen and expanding the true lumen (TL), the distal uncovered stent component enables the stabilization of the dissection flap within the aortic arch and descending aorta thereby stabilizing the structure of the aortic wall and promoting healing of the wall. The proximal PTFE felt component is used to buttress and strengthen the aortic tissue in preparation to perform the conventional polyester graft to aorta anastomosis. The entire procedure is performed utilizing an open chest approach while on cardiopulmonary bypass.

The AMDS is mounted on a delivery system as shown in Figure 2. The delivery system consists of a tip, inner rod with the loaded device, a protective sheath, a deployment mechanism and a handle component (Figure 2). The AMDS is located just proximal to the tip of the delivery device and stretched on the inner rod to reduce its profile. The braided stent component of the AMDS is then held in place by an ePTFE suture circumferentially knotted around the loaded stent to stabilize the braided stent onto the inner rod. The suture restrains the stent allowing for sequential deployment of the stent component from proximal to distal direction. The ePTFE suture is then passed through a lumen in the inner rod and connected to a green cap on the deployment mechanism as shown in Figure 2. The red cap on the delivery system is non-functional. The handle component is equipped with an end port for guide wire exit and for flushing the wire lumen with saline during the preparation of the device immediately prior to use. The most proximal aspect of the stent and the graft component is constrained by a 6 cm transparent plastic protective sheath that reduces the profile of the proximal graft component. This plastic sheath allows for atraumatic insertion of the proximal aspect of the stent into the aorta.


Figure 1- AMDS Hybrid Prosthesis (40mm shown)


Figure 2- AMDS Hybrid Prosthesis (40mm shown)

Table 1. Characteristics of the AMDS

REF	Device	Length Free State (mm)	Max Length (mm)*	Catheter Length (mm)	Introducer	Guidewire
AMDS 40	AMDS 40	95	208	527	Soft silicone tube diameter <10mm	0.035"
AMDS 40c	AMDS 40 (Classic)	95	208	527	Soft silicone tube diameter <10mm	0.035"
AMDS 4030	AMDS 40-30 Tapered	92	210	527	Soft silicone tube diameter <10mm	0.035"
AMDS 4030c	AMDS 40-30 Tapered (Classic)	92	210	527	Soft silicone tube diameter <10mm	0.035"
AMDS 55	AMDS 55	121	231	565	Soft silicone tube diameter <10mm	0.035"
AMDS 55c	AMDS 55 (Classic)	121	231	565	Soft silicone tube diameter <10mm	0.035"
AMDS 5540	AMDS 55-40 Tapered	112	225	546	Soft silicone tube diameter <10mm	0.035"
AMDS 5540c	AMDS 55-40 Tapered (Classic)	112	225	546	Soft silicone tube diameter <10mm	0.035"

* Represents the maximum stent length when the stent expansion is at the smallest aorta diameter according to the instructions for use of the particular AMDS referenced here.

2.0 Mechanism of Action

Two anatomical components, the initial entry tear and distal extension of the intimal flap contribute to the acute and chronic complications associated with acute DeBakey type I dissections. Therefore, to accomplish successful repair, the initial tear has to be excluded/ removed and the distally extending dissection flap needs to be re-expanded and supported to reestablish flow through the true lumen and promote remodeling. The AMDS uses the distal uncovered component of the device implanted distal to the aortic anastomosis, to expand the true lumen and support the intimal flap, thereby promoting remodeling in the aortic arch. In addition, the stent will expand the TL and improve the blood flow through the aorta and its tributaries. Because the stent component is uncovered, it will allow uninhibited flow to the aortic side branches. The proximal PTFE felt graft component is used to strengthen the distal aortic anastomosis created between a conventional polyester graft and the transected aorta. The proximal PTFE component serves to seal the aortic layers and fully close the false lumen at the distal anastomosis. The PTFE felt is not intended as replacement for aortic tissue, but it will strengthen the aortic tissue to create a more stable and hemostatic anastomosis between the conventional polyester graft and the aorta. The main function of the PTFE felt component is to assist in closing the false lumen at the site of the conventional graft to aorta anastomosis.

3.0 Packaging and Sterilization

Product is secured in a tray with a Tyvek® sealed lid. The product is transported inside of a single unit carton. Note: The single unit carton is not sterile. Only the AMDS may be introduced to the sterile field. The AMDS is sterilized by ethylene oxide, is supplied sterile, and must not be re-sterilized. The Tyvek® seal on the tray must be intact. Any damage to the seal renders the prosthesis non-sterile. In the event of damage to the primary packaging, the product must not be used and should be returned immediately to the supplier. The AMDS is sterilized and non-pyrogenic.

4.0 Magnetic Resonance Imaging (MRI) Safety Information

MR Conditional

Non-clinical testing and MRI simulations were performed to evaluate the entire family of AMDS. Non-clinical testing and electromagnetic simulations demonstrated that the entire family of this product is MR Conditional. A patient with an implant from this family can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 2,000-gauss/cm (20-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, the AMDS implant is expected to produce a maximum temperature rise of 5.0°C after 15-minutes of continuous scanning (i.e., per pulse

sequence). In non-clinical testing, the image artifact caused by an implant from AMDS extends approximately 12-mm from this device when imaged with a gradient echo pulse sequence and a 3-Tesla MR system. The lumen of this stent cannot be visualized on the gradient echo pulse sequence.

5.0 Intended Use/Indications for Use

The AMDS is indicated for the treatment of Acute DeBakey Type I Dissections. AMDS is intended for aortic repair, aortic remodeling and re-expansion of the intimal flap within the ascending aorta, aortic arch, and into the descending aorta for patients with acute DeBakey Type I aortic dissections undergoing open surgical repair within 0-14 days of diagnosis.

6.0 Contraindications

The AMDS should not be implanted in patients who exhibit sensitivity to PTFE or Nitinol (e.g. Nickel or Titanium), patients with mycotic aneurysms, aortic fistulous communication with non-vascular structures, pediatric patients < 18 years, or patients with uncontrolled systemic infection.

7.0 Cautions and Warnings

1. The AMDS PTFE felt should not be used as a graft to replace the aorta. The PTFE felt is solely used to buttress and strengthen the anastomosis between the polyester graft and the aorta.
2. All dissected aortic tissue should be resected and removed proximal to the AMDS. The conventional aortic graft must be anastomosed to the aorta incorporating the conventional graft, aortic wall and the PTFE felt in each suture bite.
3. Because of the discrepancy between the size of the AMDS graft component and the diameter of the various aortas treated at the site of AMDS attachment, the excess aortic tissue may create folds that would serve as entry point of blood flow into the false lumen creating a persistent false lumen. To avoid this, all aortic tissue proximal to the AMDS should be resected and all the gaps meticulously sealed and incorporated into the conventional graft to the aorta-AMDS complex anastomosis as described in this document.
4. Once the PTFE felt has been attached to the aorta, the supra-aortic branches should be visually inspected by the operator to avoid any inadvertent blockage of these vessels by the PTFE felt component.
5. Do not preclot.
6. Do not use beyond the indicated expiration date.
7. Do not re-sterilize. For single use only. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, may result in deterioration of health or death of patients. Reuse, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection. Contamination of the device may lead to injury, illness or death of the patient.
8. Clamping may damage any vascular prostheses. Do not use a clamp on the PTFE felt component, stent component or the delivery system.
9. The conventional aortic graft must be anastomosed to the aorta incorporating the conventional graft, aortic wall and the PTFE felt in each suture bite. DO NOT apply tension to the AMDS collar or proximal graft while suturing the anastomosis, as this may result in stent malposition. When completing the proximal anastomosis, leave the ascending graft long enough to avoid putting tension on the AMDS at the distal anastomosis.
10. Professional use only. The AMDS procedure shall only be performed by physicians trained in cardiac surgery and hemiarch repair.
11. Do not oversize. The braided design of the AMDS stent does not require oversizing. Some degree of oversizing is inherently built into the sizing guidance as device selection is based on total aortic diameter adventitia-to-adventitia, while the native true lumen is known to be smaller than this. Device sizing must be selected by the physician based on the patient-specific anatomy, according to the AMDS Sizing Chart in Table 2. Sizing must be based on measurement of the outer aortic diameter (adventitia to adventitia) based on preoperative CT scan imaging. Implantation of the AMDS in anatomy beyond the sizing recommendations (e.g., oversizing) may result in unexpected device conformability (such as stent narrowing).
12. Use caution during removal of the delivery system to ensure the tip of the delivery system is free from the distal stent of the AMDS. Do not forcefully pull on the delivery system after deployment of the AMDS. If there is resistance during the removal of the delivery system, the curled tip may be straightened by placing a guide wire through the device, thereafter slowly removing the device.
13. Do not deploy the AMDS in the false lumen of the aortic dissection.
14. The diameter of the transected aorta in the location of the distal graft- aorta anastomosis should be 32 mm or larger to avoid plication of the PTFE felt within the aorta.
15. Patients with subacute or chronic dissections (> 14 days after the index event).
16. Extensive thrombus or calcification in the aortic arch may increase the risk of stroke.
17. Excessive aortic tortuosity may preclude safe passage of the device within the aorta.
18. Pregnant or breastfeeding women need to be counseled on the radiation risks associated with follow-up CT scans.
19. The safety and effectiveness of the AMDS has not been clinically studied for use in the patients with Marfan Syndrome or connective tissue disorders.
20. Do not remove the protective sheath until the device is advanced in the vessel. Premature removal of the protective sheath may result in device damage or pre-activation, which may result in patient injury.
21. While placing the tacking sutures in place, use care to ensure that the tacking suture is not looped around the deployment suture as this may cause damage to the deployment suture.
22. No specific anti-platelet or anticoagulant therapy should be instituted in accordance with standard practice following surgical repair of ATAD and at the discretion of the treating medical team.

8.0 Potential Adverse Events

- Allergic reaction (to contrast media, antithrombotic agents, prosthesis material)
- Amputation
- Aortic enlargement (e.g., persisting flow in the false lumen)
- Aortic rupture
- Arterial or venous thrombosis and/or pseudoaneurysm
- Cardiac complications and subsequent problems (e.g., arrhythmia, tachycardia, tamponade, myocardial infarction, hypotension, hypertension)
- Cardiac failure (e.g., congestive heart failure)
- Death
- Dissection, perforation, or rupture of the aortic vessel & surrounding vasculature
- Embolism (e.g., thromboembolism)
- Fistula (e.g., aorto-esophageal, aorto-tracheal, aorto-bronchial)
- Gastrointestinal complications and subsequent problems (e.g., visceral ischemia/infarction, nausea, vomiting)
- Hemorrhage/bleeding
- Hepatic failure
- Infection (e.g., local, systemic, prosthesis) or fever
- Ischemia or infarction (e.g., cerebral, visceral, renal, organ, peripheral)
- Neurological complications and subsequent problems (e.g., transient ischemic attack (TIA), stroke, neuropathy)
- Occlusion (venous or arterial, incl. prosthesis occlusion)
- Pain and inflammation
- Pulmonary complication (e.g., edema, embolism, pneumonia, respiratory failure)
- Renal insufficiency
- Spinal cord ischemia (SCI) incl. paraparesis and paraplegia
- Stenosis (arterial or venous)
- Wound complications and subsequent problems (e.g., dehiscence, infection)

Sepsis

9.0 Device sizing

The AMDS is supplied sterile in 8 configurations. While the diameter of the proximal graft component in both devices is identical, the stent diameters supplied are 40mm straight, 40mm tapered (30mm diameter distal end), 55mm straight and 55mm tapered (40mm diameter distal end). The AMDS has the capability to change its diameter as a function of its length, i.e. the device lengthens or shortens as its diameter is reduced or increased. This allows for treatment of patients with aortic diameters between 20-45mm according to the table below. The aorta is sized based on a contrast enhanced CT scan or MR- Angiography of the chest, measuring a) the largest diameter of the aortic arch as defined between the distal border of the innominate trunk and the proximal descending thoracic aorta, b) the diameter of mid-descending aorta at the levels of T6, and c) the length of the aorta (aortic length*) as defined by the distance between the planned proximal location of implantation of the AMDS to the aorta immediately proximal to the Celiac trunk. For measurements, two and/ or three-dimensional images may be used. The aortic diameter is defined by adventitia-to-adventitia (outer to outer) measurements. To ensure that the AMDS is not placed across the abdominal visceral vessels, prior to the implantation of the AMDS, the operator has to ensure that the aortic length measures to be equal or longer than that stated for the unit lengths at minimum IFU diameter according to table 2. For example, if the measured proximal aortic diameter is 25mm and the distal aortic diameter is 20mm, the measured aortic length* should be no shorter than 195mm.

Table 2. Device sizing table.


REF	Device	Proximal Aorta Diameter (mm)	Distal Aorta Diameter (mm)	Unit Length at Minimum IFU Diameter (mm)	Unit Length at Maximum IFU Diameter (mm)
AMDS 40	AMDS 40	20-35	25-35	208 mm	155 mm
AMDS 40c	AMDS 40 (Classic)	20-35	25-35	208 mm	155 mm
AMDS 4030	AMDS 40-30 Tapered	20-35	20-24	210 mm	170 mm
AMDS 4030c	AMDS 40-30 Tapered (Classic)	20-35	20-24	210 mm	170 mm
AMDS 55	AMDS 55	36-45	36-45	231 mm	195 mm
AMDS 55c	AMDS 55 (Classic)	36-45	36-45	231 mm	195 mm
AMDS 5540	AMDS 55-40 Tapered	36-45	27-35	225 mm	190 mm
AMDS 5540c	AMDS 55-40 Tapered (Classic)	36-45	27-35	225 mm	190 mm

* Aortic length= the distance between the planned implantation site of the proximal component of the AMDS to the aorta immediately proximal to the Celiac trunk
Classic codes are not approved or marketed in Australia.

10. Implantation

10.1 Preparation of the Aorta and Distal Anastomosis

The ascending aorta (including the aortic arch if aneurysmal) is transected and removed in a routine standard of care manner utilizing hypothermic circulatory arrest. The operator must ensure to leave at least 10 mm (1.0 cm) aortic tissue proximal to the native arch branch vessels that will be reserved in order to avoid inadvertent coverage of supra-aortic branches with the covered component of the AMDS. Initially the lesser curvature should be left longer to simplify the attachment of the AMDS to the aorta. The excess aortic tissue should be resected up to the level of the AMDS PTFE graft after the AMDS is secured to the aorta and prior to performing the anastomosis to the conventional graft. The diameter of the proximal anastomosis should be tailored to measure between 30-40 mm in order to avoid plication of the PTFE felt within the aorta.

10.2 Preparation of the AMDS


Visually inspect the device packaging. If the package has been previously opened or the seal is damaged, discard and replace with a new package. The entire device and delivery system must be pre-soaked and flushed with saline to ensure adequate lubrication and should not be allowed to dry out after soaking.

10.3 Use of a GuideWire

The AMDS is designed to be delivered over a stiff 0.035" guide wire such as a Lunderquist guide wire, if desired. The entry point for the guide wire is at the tip of the delivery system. The wire will exit the delivery system in the handle component. This allows the system to track into position along the natural curvature of the aorta. Ideally the guidewire should be inserted using a transfemoral technique prior to placing the patient on cardiopulmonary bypass or while on cardiopulmonary bypass prior to inducing hypothermic arrest. The position of the guidewire should be confirmed within the true lumen using TEE or IVUS.

10.4 Positioning of the Device


More than 1.0 cm of aortic tissue should be left in place proximal to the innominate artery or any of the arch branches that are intended to be preserved. The AMDS is delivered into the true lumen through the open distal aorta. This can be done over a guide wire to ensure access to the correct lumen (true lumen) in the dissected aorta. The AMDS is gently inserted into the aorta until the PTFE felt collar of the device is in the same plane as the transected distal aorta. The device is symmetrical and therefore does not require rotational orientation.

10.5 Deployment Sequence


After the positioning of the AMDS in the aorta has been completed and the edge of the felt collar is at the level of the transected edge of the aorta, the protective sheath is retracted to fully expose the PTFE felt. The protective sheath is removed and discarded. The PTFE felt is stabilized by placing interrupted or running polypropylene sutures, attaching the PTFE felt collar to the proximal-most inner aspect of the transected aorta. Use of a PTFE felt strip external to the aorta is highly recommended to avoid excessive tension on or tearing of the aorta. Once the PTFE felt has been attached to the aorta, the supra-aortic branches should be visually inspected by the operator to avoid any inadvertent blockage of these vessels by the PTFE felt component. The inspection is best performed by using a sterile scope.

While an assistant stabilizes the PTFE felt to the aorta using gentle vascular pickups, the operator unscrews the green cap in a counterclockwise fashion releasing the green cap. The cap and the ePTFE wrapping suture can now be pulled back to release the braided stent and allow it to expand. The braided stent opens from a proximal to a distal direction on the delivery system. During the release of the wrapping suture expanding the stent, the operator needs to apply firm but gentle traction to the delivery system. Once the stent is released and disengaged from the delivery mechanism, the delivery system is removed making sure the tip of the delivery system is free from the distal end of the stent. If the operator experiences any resistance in removing the delivery system after the AMDS is deployed, a guidewire may be placed through the delivery system to uncurl the tip prior to removing the delivery system.

10.6 Anastomoses

Once the delivery system is removed, we recommend preparing the aortic anastomotic site by initially completely securing the aortic layers and sealing the false lumen by running a continuous suture layer circumferentially around the transected aorta, attaching the inner PTFE felt (AMDS) through the aorta and to the outer PTFE felt. There should be no folding of the inner felt as it may allow blood flow to find its way to the false lumen once flow is established. Once this suture line is complete, any of the approved aortic replacement polyester grafts may be chosen and the graft to aorta distal anastomosis performed in a conventional way, assuring the AMDS graft collar is used to buttress the polyester graft to aorta anastomosis from within the proximal-most aspect of the transected aorta. All dissected aortic tissue should be removed proximal to the AMDS. The use of a PTFE felt strip external to the anastomosis is highly recommended as described above. The anastomosis between the conventional polyester graft and the AMDS-aorta complex should be performed with each suture bite incorporating the conventional graft and the PTFE felt- aorta- PTFE felt (AMDS) complex creating a meticulously tight seal between the PTFE felts and the dissected aorta. This technique in essence allows the operator to anastomose the conventional graft to the AMDS “washer” created from external felt- aorta- internal felt (AMDS) to create a sealed and secure anastomosis.

11.0 Storage

The AMDS, in its protective packaging, should be stored at room temperature, specifically not less than 0°C and not more than 25°C, and in a dry location.

12.0 Disposal

At the end of the procedure care must be taken to ensure safe disposal of the AMDS delivery system. The surgical team must ensure local and national regulatory requirements for the disposal of contaminated clinical waste products is adhered to.

13.0 References

Tyvek® Du Pont Registered Trademark. *AMDS is a trademark of Artivion, Inc. or its subsidiaries

14.0 Symbols used on labeling

	Manufacturer		Use By Date		Keep Dry		Patient Name or Patient ID
	Date of Manufacture		Consult Instruction for Use		Temperature Limit		Date of implantation
	Catalog Number		Sterilized Using Ethylene Oxide		MR Conditional		Name and Address of the implanting healthcare institution/provider
	Batch Code		Do not reuse		Information website for patients		
	Medical Device		Do Not Use If Package Is Damaged		UDI as AIDC format		CE Mark

AMDS Hybrid Prosthesis Implant Card

UDI

A patient with an AMDS implant can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 2,000-gauss/cm (20-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode


Under the scan conditions defined, the AMDS implant is expected to produce a maximum temperature rise of 5.0°C after 15-minutes of continuous scanning (i.e., per pulse sequence). In non-clinical testing, the image artifact caused by AMDS extends approximately 12-mm from this device when imaged with a gradient echo pulse sequence and a 3-Tesla MR system. The lumen of this stent cannot be visualized on the gradient echo pulse sequence.

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Tyvek® Du Pont Tescilli Ticari Marka. AMDS, Artivion, Inc. veya bağlı şirketlerinin ticari markasıdır.