CAUTION – Investigational device. Limited by Federal (or United States) law to investigational use.
Exclusively for clinical investigation.

For eIFU, visit: www.artivion.com/eifu/arcevo

A printed IFU will be provided within seven days upon request through any of the
contact means listed below:

Phone: +1 888.427.9654 • Fax: 770.590.3753 • Email: customerservice.us@artivion.com

Artivion, Inc.
1655 Roberts Blvd. NW
Kennesaw, GA 30144 USA
+1 (888) 427-9654
www.artivion.com

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1.0 OVERVIEW OF THE ARCEVO™ LSA HYBRID STENT GRAFT SYSTEM

1.1 INTENDED USE

The Arcevo™ LSA Hybrid Stent Graft System is intended for the surgical treatment of patients with dissection and/or aneurysm involving the aortic arch and the descending thoracic aorta, with or without involvement of the ascending aorta. The Arcevo™ LSA Hybrid Stent Graft System is intended for attachment to a proximal surgical graft (not supplied). Patients who need additional stent graft coverage distally, as determined by a physician, may require extension with the Medtronic® Valiant™ thoracic stent graft with the Captivia™ delivery system (also referred to as ‘Valiant™ Captivia™’; not supplied) after the initial index procedure using the Instructions for Use provided for that device. The Arcevo™ LSA Hybrid Stent Graft System is intended to be used by trained cardiothoracic surgeons with experience in performing aortic arch procedures in a surgical suite or in the operating room.

1.2 DEVICE DESCRIPTION

The Arcevo™ LSA Hybrid Stent Graft System (also referred to as ‘Arcevo™ LSA’) shown in **Figure 1** consists of an implantable aortic stent graft with left subclavian artery (LSA) branch which is preloaded onto a delivery system. The Arcevo™ LSA is designed to be introduced into the transected aorta (with or without the use of an aortic guidewire), and it is strongly recommended to position the LSA branch into the LSA vessel using a guidewire. The implant is crimped onto the delivery system using two polyethylene terephthalate (PET) Wraps, which are sutured to the Arcevo™ implant and laced with the Aortic and LSA Deployment Wires to constrain the implant during insertion and positioning within the aorta. The Arcevo™ LSA implant is deployed by first retracting the Aortic Deployment Wire, then the LSA Deployment Wire, to unlace each Wrap. The Wraps remain secured between the stent graft and the aortic wall post-deployment. The LSA Deployment Wire is retracted into the handle prior to delivery system withdrawal. The implant collar is then sutured to a surgical graft (not supplied) placed in the aortic arch to complete the total arch repair procedure. The Arcevo™ LSA System is packaged in a double-sterile barrier consisting of an inner and outer pouch secured inside a shelf box and provided sterile (via ethylene oxide) for single-use only. The Arcevo™ LSA is available in multiple configurations as shown in **Table 1**. The aortic stent graft is available in a variety of diameters and either 120 mm or 130 mm lengths, with an LSA branch available in either an 11 mm or 15 mm diameter, and a 40 mm length.


Figure 1. Arcevo™ LSA Hybrid Stent Graft System


Figure 2. Arcevo™ LSA Implant

Table 1. Arcevo™ LSA Device Configurations – see Figures 2 and 3 for Implant Image and Dimensions

Arcevo™ LSA Model	[A] Aortic Diameter (mm)	[B] LSA Diameter (mm)	[C] Nominal Aortic Length* (mm)	[D] Nominal LSA Length* (mm)	[E] Collar Diameter (mm)	[F] LSA Branch Angle (°)	[G] Collar to Branch Length (mm)
ARC-A24-B11	24	11	120	40	50	70	10
ARC-A26-B11	26	11					10
ARC-A28-B11	28	11					12
ARC-A28-B15	28	15			55		10
ARC-A30-B11	30	11					12
ARC-A30-B15	30	15					10
ARC-A33-B11	33	11	130	40	60	12	
ARC-A33-B15	33	15				10	
ARC-A36-B11	36	11				12	
ARC-A36-B15	36	15				10	
ARC-A40-B15	40	15				10	

*Nominal aortic lengths range 119-124 mm for 120 mm devices and 132-134 mm for 130 mm devices. Nominal LSA branch lengths range 40-44 mm.

The Arcevo™ LSA implant (**Figure 2**) is designed with an aortic stent graft body and a stented LSA branch. The implant is made of self-expanding nitinol aortic stents attached to polyethylene terephthalate (PET/polyester) graft material with ultra-high molecular weight polyethylene (UHMWPE) sutures. The proximal end of the stent graft has a conical shaped collar which is approximately 15 mm in length for anastomosis to the native aorta and a proximal surgical graft (not supplied). The aortic body includes five platinum-iridium radiopaque markers: three at the proximal end of the stent graft around the base of the LSA branch, one along the aortic graft body to indicate an overlap distance of 50 mm, and one at the distal end. The LSA branch also has one platinum iridium marker at the distal tip of the branch.


Figure 3. Arcevo™ LSA Implant Dimensions

The Arcevo™ LSA Delivery System (Figure 4) includes aortic and LSA sections which are connected to the handle for insertion, positioning, and deployment of the Arcevo™ LSA implant into the aorta. Each section of the delivery system is comprised of polyether ether ketone (PEEK) guidewire tube ending in an atraumatic acetal polyoxymethylene (POM) copolymer tip designed to guide the delivery system into the thoracic aorta and LSA vessel. The aortic guidewire tube is threaded with twenty-four polybutylene terephthalate (PBT) pearls which provide a surface for crimping of the aortic stent graft. The LSA guidewire tube includes a single pearl for LSA branch securement. A flexible stainless steel Shaping Wire is provided loaded within the delivery system to add delivery system support during implant pre-shaping, insertion, and deployment if an aortic guidewire is not used. As mentioned previously, deployment occurs by retracting the Aortic Deployment Wire from the Wrap constraining the aortic stent graft by pulling on the trigger and end of handle. Subsequently, the LSA branch is deployed by pulling the LSA Tab to retract the LSA Deployment Wire from the LSA Wrap constraining the LSA branch.


Figure 4. Arcevo™ LSA Delivery System

2.0 INDICATIONS FOR USE

The Arcevo™ LSA Hybrid Stent Graft System is indicated for use in patients with dissection and/or aneurysm involving the aortic arch and the descending thoracic aorta, with or without involvement of the ascending aorta.

3.0 CONTRAINDICATIONS

The Arcevo™ LSA implant is contraindicated in patients with known sensitivities or allergies to nitinol, platinum-iridium, polyester, and/or polyethylene, in patients with systemic infection, and in patients with endocarditis or active infection of the aorta.

4.0 WARNINGS

- **Professional use only.** The Arcevo™ LSA procedure shall only be performed by physicians trained in cardiac surgery and total arch replacement. Patient selection and device sizing must be chosen by the physician based on the patient-specific anatomy using the recommended sizing guidelines. Sizing must be based on measurement of the aortic diameter based on preoperative imaging. Implantation of the Arcevo™ LSA in anatomy beyond the recommended guidelines may result in patient injury, treatment failure requiring additional intervention, and possible death.
- **For single use only. DO NOT reuse, reprocess, or re-sterilize.** Reuse, reprocessing, or re-sterilization may compromise structural integrity of the device and/or lead to device failure which may result in deterioration of health or death of patients. Reuse, reprocessing, or re-sterilization may lead to contamination of the device and cause patient infection or transmission of infectious disease leading to injury, illness, or death of patient or user.
- **DO NOT** use the device if the labeled expiration date has passed, the packaging is opened or damaged, and/or the product has been dropped, damaged, or otherwise inappropriately handled. The sterile inner pouch must be intact at the time of use; damage to the seal renders the Arcevo™ LSA non-sterile. If the device or packaging is damaged, return the product to the manufacturer.

5.0 PRECAUTIONS

- Safety and effectiveness have not yet been established for the Arcevo™ LSA device. Precaution should be taken in the following patient populations:
 - Patients with history of bleeding diathesis or coagulopathy, who may have an increased risk of severe bleeding or hematoma.
 - Patients with severe angulation, circumferential thrombus, or severe calcification in the distal landing zone may increase the risk of migration or endoleak formation. To achieve distal seal in an aneurysm, a distal aortic neck length of at least two distal springs of the stent graft is required.
- In order to enhance organ protection during surgery using the frozen elephant trunk technique, carefully consider individualization of treatment in the following areas: hypothermic circulatory arrest and degrees of hypothermia, cerebral perfusion (antegrade or retrograde), and distal body perfusion.
- Adherence to the device sizing guidelines is recommended – tortuosity, landing zone, and thrombus, and/or calcium at the implantation sites may lead to endoleak, migration at the seal zone(s), and/or device damage.

- Vessel over-distension and damage may be caused by excessive oversizing of the stent graft in relation to the diameter of the vessel.
- Patients with aortic dissection at the origin of the visceral and renal arteries should be assessed preoperatively; if vessel perfusion depends on a patent false lumen, treatment using the FET technique may cause renal, visceral, and/or peripheral ischemia or infarction.
- Severe atherosclerosis at the distal landing zone may result in increased risk of spinal cord ischemia as a consequence of thromboembolism of intercostal arteries or branch vessels.
- Severe calcification of the native aortic tissue at the site of anastomosis may lead to bleeding or difficulty in suturing.
- Patients should be anticoagulated per institutional protocol. Insufficient anticoagulation during the implantation procedure may result in thrombosis and embolism. Leakages through the graft fabric (endoleak type IV) which can be visualized during or immediately after the implantation may be caused by intraoperative anticoagulation. Any endoleak left untreated during the implantation procedure must be carefully monitored after implantation.
- Clamping may damage the prosthesis. Atraumatic clamps, ideally with soft shod jaws, should be used with a minimum application of force. Excessive force should be avoided as it will damage the polyester fibers or nitinol stents.
- Incomplete or inaccurate anastomosis may result in bleeding and/or the need for re-exploration.
- Care should be taken to include all dissected membranes of the aortic wall in the anastomosis. In case of very arteriosclerotic and friable tissue, the use of external felt to support the anastomosis may minimize the risk of bleeding.
- Excessive tension on the anastomosis and the graft material should be avoided to prevent strain on the suture line.
- If de-airing is required, the smallest needle possible should be used. DO NOT use cutting needles that may result in blood leakage.

6.0 POTENTIAL ADVERSE EVENTS AND UNDESIRABLE SIDE EFFECTS

Potential risks that may be associated with the total arch replacement (TAR) procedure include, but are not limited to:

- Anemia requiring ≥ 2 units of blood product in a transfusion
- Atrial fibrillation
- Biologic response (e.g., antithrombotic, systemic inflammatory response syndrome [SIRS])
- Cardiac failure (e.g., aortic regurgitation, arrhythmia, congestive heart failure, myocardial infarction, pericardial effusion, tamponade)
- Death
- Fistula (e.g., aorto-esophageal, aorto-bronchial)
- Hematoma
- Hemorrhage/bleeding
- Hypotension/hypertension
- Infection (e.g., local, systemic, prosthesis) or fever
- Infectious disease transmission
- Inflammation (e.g., edema, pancreatitis, gout, etc.)
- Internal bleeding
- Ischemia-reperfusion injury
- Recurrent laryngeal or phrenic nerve palsy
- Liver failure
- Multi organ dysfunction syndrome or failure
- Gastrointestinal system complications and subsequent problems (e.g., visceral ischemia/infarction, nausea and vomiting)
- Neurological complications (e.g., edema, headache, transient ischemic attack (TIA), change in mental status, paralysis, spinal cord ischemia (SCI) including paraparesis and paraplegia, seizure, spinal cord infarction, stroke) or temporary effects on and/or damage to neural tissue
- Pleural effusion requiring drainage (percutaneous thoracentesis or tube thoracostomy)
- Respiratory failure requiring tracheostomy or prolonged ventilation
- Pulmonary complication (e.g., edema, embolism, interstitial syndrome, mediastinitis, pneumonia, pneumothorax, respiratory depression)
- Renal insufficiency and/or renal failure
- Sepsis
- Shock
- Vascular trauma
- Vasoplegic syndrome
- Wound healing complications (e.g., dehiscence, bleeding, infection)

Potential risks that may be associated with Arcevo™ LSA device or implantation procedure include, but are not limited to:

- Allergic reaction to prosthesis material(s)
- Amputation due to upper limb ischemia
- Aortic enlargement (e.g., aneurysm, false lumen)
- Aortic aneurysm rupture
- Dissection/perforation/rupture of the aorta, LSA, and/or vasculature, distal stent induced new entry (dSINE) tears
- Embolism (e.g., thromboembolism, microembolism)
- Endoleak, distal anastomotic new entry tear (DANE) or distal anastomotic leak
- Ischemia or infarction (e.g., cerebral, visceral, renal, organ, peripheral)
- Occlusion (venous or arterial, including prosthesis occlusion)
- Post-implant syndrome
- Prosthesis material failure (including breakage of metal portion of device, collapse, extrusion, erosion, and/or prosthesis rupture)
- Pseudoaneurysm
- Stenosis (arterial or venous including prosthesis stenosis)
- Stent graft migration (retrograde in case of aneurysm)
- Stent graft infection
- Thrombosis (including prosthesis thrombosis)

7.0 MAGNETIC RESONANCE IMAGING SAFETY AND COMPATIBILITY INFORMATION

A patient with the Arcevo™ LSA (with or without a Medtronic Valiant™ thoracic stent graft) can be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient. If information about a specific parameter is not included, there are no conditions associated with that parameter.

	Static Magnetic Field Strength (B_0)	1.5T or 3T
	Maximum Spatial Field Gradient	30-T/m (3,000-gauss/cm)
	RF Excitation	Circularly Polarized (CP)
	Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
	Maximum Head SAR	3.2 W/kg (Normal Operating Mode)
	Limits on Scan Duration	2.0 W/kg whole body average SAR for 60 minutes of continuous RF (a sequence or back-to-back series/scan without breaks)
	MR Image Artifact	The implant(s) may produce an image artifact.
	MR Conditional	

Patients who have other MR Conditional devices can be scanned as long as all the MR Conditional scan parameters for each of the devices are met. DO NOT conduct an MRI scan if any conditions for safe scanning for any device cannot be met.

8.0 PATIENT SELECTION GUIDELINES AND DEVICE SIZING

WARNING: DO NOT use the Arcevo™ LSA in patients with unsuitable LSA anatomy or disease states, including aneurysmal, or stenotic LSA, extreme LSA take-off angle, anomalous LSA anatomy, and/or low left vertebral artery take-off location.

8.1 PATIENT SELECTION

Before use of the Arcevo™ LSA Hybrid Stent Graft System, physicians should carefully consider the benefits and risks for each patient with respect to patient age, life expectancy, comorbidities, ability to tolerate general anesthesia, and suitability for open surgical repair. The Arcevo™ LSA is not recommended in patients who are unable to undergo or not compliant with the necessary pre- and post-operative imaging. Patients who cannot tolerate contrast agents necessary for pre- and post-operative follow-up imaging are not recommendable. In case of acute dissection, clinical judgement must be made per specific patient and clinical circumstance.

The eligibility criteria for the ARTIZEN study are provided below:

Inclusion Criteria:

General Inclusion Criteria

1. ≥18 years of age or ≤80 years of age (male or female) at time of surgery
2. Patient has one of the following indications for open surgery based on CTA completed within 90-days of informed consent:
 - Acute, subacute, or chronic dissection that involves the aortic arch and the descending thoracic aorta, with or without involvement of the ascending aorta
 - Aneurysm that involves the aortic arch and the descending thoracic aorta, with or without involvement of the ascending aorta
3. Patient, or patient's legally authorized representative (in the secondary arm only), provides written informed consent prior to any study procedures
4. Patient's surgery occurs within 90 days of informed consent

Anatomical Inclusion Criteria

5. Aortic diameter at the intended Arcevo™ LSA anastomosis site is ≥ 21 mm
6. For Aneurysm with distal sealing (i.e., single stage procedure), aortic diameter at the intended distal sealing zone is between 20-36 mm
7. For Dissection with distal sealing, aortic diameter at the intended distal sealing zone is between 22-40 mm
8. LSA branch does not require additional stenting further into the LSA (beyond the Arcevo™ LSA)
9. The intended LSA sealing zone has a diameter between 8.5 mm and 14.0 mm with a length of at least 10 mm
10. Absence of dissection, aneurysm, or stenosis in the intended LSA sealing zone
11. The intended LSA branch position does not interrupt flow to any branch vessel (e.g., left vertebral artery)
12. Patient does not have CT evidence of extreme arch or LSA angulation precluding safe passage of the Arcevo™ LSA and delivery system
13. LSA take-off angle between 15° and 90°
14. For patients requiring planned extension, there is a ≥ 2 cm distal seal zone for the TEVAR device without a severely angulated descending aorta
15. For patients requiring planned extension, patient has suitable iliac artery anatomy for safe passage of the TEVAR delivery system

Exclusion Criteria:

General Exclusion Criteria

1. Patient is pregnant, or planning to become pregnant during the course of the study; individuals of child-bearing potential must agree to use acceptable methods of contraception during the study
2. Patient has another medical condition (aside from the arch disease) that, in the opinion of the investigator, reduces the patient's life expectancy to < 2 years
3. Patient has an existing aortic stent graft device in the descending aorta that would interact with Arcevo™ LSA
4. Patient has a medical, social, or psychological problem that, in the opinion of the investigator, could impede the patient's ability to return for follow-up
5. Patient is unwilling or unable to comply with the follow-up schedule
6. Patient is institutionalized due to administrative or judicial order
7. Patient is unwilling to accept blood transfusion or blood product
8. Patient is currently participating in another interventional clinical study which includes treatment with another investigational product (e.g., device, pharmaceutical or biologic)

Medical Exclusion Criteria

9. Patient is unfit for open surgical repair involving circulatory arrest
10. Patient is in extreme hemodynamic compromise requiring cardiopulmonary resuscitation or substantial inotropic support prior to surgery
11. Patient has an active systemic infection
12. Patient has endocarditis or active infection of the aorta
13. Patient has a freely ruptured aorta
14. Patient has a history of a bleeding disorder (e.g., hemophilia)
15. Patient has current end-stage renal disease (e.g., GFR <30mL/min)
16. Patient has a known allergy to Arcevo™ LSA components and TEVAR device components (if required) (nitinol, polyester, platinum-iridium, or polyethylene)
17. Patient has uncontrollable anaphylaxis to iodinated contrast or other inability to obtain CT angiograms during follow-up
18. Patient has acute coronary malperfusion
19. Patient has symptomatic visceral malperfusion

8.2 DEVICE SIZING AND PLACEMENT

The Arcevo™ LSA implant is designed to treat patients with aortic aneurysms (aortic diameters 20-36 mm) and aortic dissections (aortic diameters 22-40 mm), with LSA diameters between 8.5-14.0 mm. The specific stent graft diameter should be selected for individualization of treatment according to the treated disease using the recommended sizing guidelines and the provided anatomical and clinical criteria. Sizing of the Arcevo™ LSA implant should be based on the complete treatment strategy, including the size of the distal landing zone of the thoracic extension stent graft device, if needed.

WARNING: Failure to size the Arcevo™ LSA implant within the recommended guidelines based on aortic disease and patient anatomy may result in patient injury; excessive oversizing may lead to aortic injury or rupture. Undersizing the Arcevo™ LSA implant or thoracic extension device may lead to device stenosis, migration and/or endoleak, or possible reintervention. Clinical judgement shall be used in determining Arcevo™ LSA stent graft oversizing to the native aortic anatomy.

The treating physician is responsible for selecting the appropriate device sizes. **Table 2** and **Figure 5** are considered to be tools for selecting the appropriate size. If a distal thoracic extension is planned, the distal end of the Arcevo™ LSA stent graft represents an adequate landing zone. The overlapping should be at least 50 mm. The chain shaped overlap marker indicates the required overlapping. Consult the Medtronic Valiant™ Thoracic Stent Graft with the Captivia™ Delivery System Instructions for Use for distal thoracic extension sizing guidelines.

Table 2. Arcevo™ LSA Device Sizing Guidelines

Arcevo™ LSA Model	All Indications:		If Distal Sealing is Intended:			
			Aortic diameter at the intended Arcevo™ LSA distal sealing zone (mm)		Minimum length of seal zone (mm)	
	Aortic diameter at the intended Arcevo™ LSA proximal anastomosis site (mm)	LSA diameter range (mm)	Aneurysm*	Dissection**		
ARC-A24-B11	≥ 21	8.5 – 14.0	20 – 22	22 – 24	25	
ARC-A26-B11	≥ 23		22 – 24	24 – 26		
ARC-A28-B11	≥ 25		24 – 25	26 – 28		28
ARC-A28-B15			25 – 28	28 – 30		
ARC-A30-B11	≥ 26		28 – 30	30 – 33	33	
ARC-A30-B15			30 – 34	33 – 36		
ARC-A33-B11	≥ 29		34 – 36	36 – 40		
ARC-A33-B15						
ARC-A36-B11	≥ 32					
ARC-A36-B15						
ARC-A40-B15	≥ 35					

*Sizing for aortic aneurysms with distal sealing were developed with approximate oversizing of 5-20% based on native aortic diameter ranges.

**Sizing for aortic dissections with distal sealing were developed with approximate oversizing up to 10% based on native aortic diameter ranges.

Contrast-enhanced computed tomographic angiography (CTA) or magnetic resonance angiography (MRA) is required to accurately assess patient anatomy prior to treatment. Measure the length of the stent graft's location at the outer aortic curvature and consider extent of inner curvature angulation to avoid the risk of stent graft kinking.

Physicians should evaluate each patient to determine if the Arcevo™ LSA implant would be appropriate for the treatment according to the inclusion and exclusion criteria above and the following considerations:

- Oversizing of the LSA branch should not exceed 30% of the native LSA vessel diameter.
- The Arcevo™ LSA collar diameter ranges from 50-65 mm; when sizing the Arcevo™ LSA implant, ensure native aortic diameter at anastomosis site is compatible with collar diameter for secure anastomosis and ideal LSA branch position. DO NOT undersize the Arcevo™ LSA implant beyond what the suture collar can accommodate at the anastomosis site.
- Extensive exclusion of the intercostal arteries by the stent graft may lead to spinal cord ischemia. Coverage of the descending thoracic aorta distal to T8 should be avoided. Staged procedures with placement of a distal thoracic stent graft should be considered in patients requiring extended coverage.


Figure 5. Acceptable Aortic Morphology for Arcevo™ LSA Patient Selection

9.0 DEVICE STORAGE, RETURNS, AND DISPOSAL

The Arcevo™ LSA Hybrid Stent Graft System, in its protective packaging, should be stored at room temperature (+5°C (41°F) to +25°C (77°F)) and kept dry. At the end of the procedure, care must be taken to ensure safe disposal of the Arcevo™ LSA delivery system. The surgical team should dispose of the Arcevo™ LSA delivery system in a biohazard waste receptacle (e.g., bag) per hospital procedure. Prior authorization from Artivion is required for the return of any product. For any questions regarding Arcevo™ LSA Hybrid Stent Graft System or for return authorization, please contact the Artivion Clinical Team at ARTIZEN@artivion.com.

10.0 EQUIPMENT

- Surgical instruments for total arch repair (scalpel, atraumatic clamps, arterial puncture needles, forceps)
- Conventional ascending aorta polyester graft
- Sterile drape
- Sterile surgical sutures for graft anastomosis to Arcevo™ LSA collar (size 3-0, 4-0, and 5-0 are acceptable)
- Strongly recommended: 0.035" (0.89 mm) guidewire intended for cardiac procedures (80 cm minimum length) for LSA branch
- Optional: 0.035" (0.89 mm) guidewire intended for cardiac procedures (150 cm minimum length) for aortic body

11.0 DIRECTIONS FOR USE

CAUTION: Preoperative planning for patient selection, device sizing, access, and placement must be performed prior to preparing the patient or unpacking the Arcevo™ LSA device. Ensure adequate access compatible with the delivery system, guidewire(s), and ancillary devices before proceeding.

11.1 PATIENT PREPARATION

11.1.1 Prepare the patient for aortic surgery per institutional protocol (e.g., cardiopulmonary bypass, hypothermic circulatory arrest, etc.).

11.1.2 Expose the access site using standard surgical technique and perform a zone 2 aortic perpendicular transection exposing the aortic arch, leaving at least 10-12 mm of aorta proximal to the LSA (**Figure 6**).

CAUTION: Failure to leave an adequate amount (10-12 mm) of aorta proximal to the LSA may result in misalignment of the Arcevo™ LSA collar with the native aorta and difficult or compromised anastomosis.


Figure 6. Zone 2 Distal Ascending Aortic Transection for Arcevo™ LSA Implantation

11.2 DEVICE PREPARATION

11.2.1 Ensure the device configuration matches the sizing determined during preoperative planning. Inspect the shelf box prior to opening. If any gross damage is observed on the box, discard the damaged unit, and obtain a new Arcevo™ LSA device. Verify that the "Use By" date has not passed.

WARNING: DO NOT use the device if the shelf box is grossly damaged. DO NOT use the device if the "Use By" date has passed.

11.2.2 OUTSIDE STERILE FIELD: Unpack the Arcevo™ LSA device by opening the shelf box and removing the outer pouch. Inspect the outer pouch for any signs of puncture, seal breach, or damage. Open the outer pouch to expose the inner pouch. Inspect the inner pouch for any signs of puncture, seal breach, or damage that may affect device sterility.

WARNING: DO NOT use the device if the outer pouch or inner pouch is grossly damaged or if the device appears displaced or damaged within its protective tray.

11.2.3 Verify that the inner pouch contains a green EO indicator label which provides confirmation that the product was EO sterilized.

CAUTION: DO NOT use the device if the EO indicator label is missing or the label color is dark purple; the Arcevo™ LSA device may not be sterile. Replace the unit with a new Arcevo™ LSA device and repeat the device preparation steps.

11.2.4 ASEPTIC TRANSFER INTO STERILE FIELD: Open the inner pouch, then aseptically transfer the device within the protective tray (with lid) to the preparation table within the sterile field.

11.2.5 INSIDE STERILE FIELD: Remove the tray lid to access the device. Remove the device from the tray by lifting at the center of the handle, taking care not to manipulate the device during transfer from the tray.

CAUTION: Manipulating the device may adversely impact the release mechanisms and lead to premature deployment or other device defects.

11.2.6 If an aortic guidewire is not to be used, pre-shape the delivery system (with the Shaping Wire still in place) to better fit the curvature of the aorta as shown in **Figure 7**.

Figure 7. Pre-shaping of aortic stent graft section with Shaping Wire in place

CAUTION: Keep both hands on the implant when pre-shaping to avoid delivery system kinking. DO NOT pre-shape the implant without the Shaping Wire. Failure to pre-shape the implant may lead to difficulty tracking the implant during insertion and positioning or damage the implant and/or delivery system.

11.2.7 If the delivery system is not pre-shaped, an aortic 0.035" (0.89 mm) guidewire is recommended for introducing the Arcevo™ LSA system into the thoracic aorta. Remove the Shaping Wire and insert the aortic guidewire by preloading the wire into the delivery system through the guidewire access. Alternatively, the guidewire can be inserted via femoral approach and fed through the Arcevo™ LSA delivery system through the Aortic Tip until it exits the handle (**Figure 8**).


Figure 8. Insertion of aortic guidewire into the Arcevo™ LSA Aortic Tip after removal of Shaping Wire

11.2.8 The Arcevo™ LSA Hybrid Stent Graft System is strongly recommended to be introduced into the LSA using a guidewire. Insert the 0.035" (0.89 mm) guidewire through the LSA guidewire lumen located at the base of the LSA tab until it exits the LSA Tip. Alternatively, if retrograde axillary arm access is used, place the guidewire through the LSA until it exits the transected aorta, then feed the guidewire through the LSA Tip until it exits the LSA guidewire lumen at the LSA Pull Tab. (**Figure 9**).

Figure 9. Inserting an LSA guidewire through the access lumen of Arcevo™ LSA delivery branch

11.3 INSERTION AND POSITIONING

11.3.1 Align the LSA branch with the location of the LSA vessel prior to inserting the Arcevo™ LSA implant into the transected thoracic aorta (**Figure 10**). The Arcevo™ LSA delivery system is not designed to be rotated once inserted into the aorta.


Figure 10. Introduce the Arcevo™ LSA device to the target position

11.3.2 Advance the Arcevo™ LSA device (with the LSA branch aligned for positioning into the LSA vessel) into the transected aorta.

CAUTION: DO NOT rotate the Arcevo™ LSA device once the device has been inserted into the aorta. Failure to align the LSA branch with the LSA vessel prior to insertion into the aorta may require removal and re-insertion of the device and may result in vessel trauma or premature deployment of the implant.

11.3.3 Advance the stented LSA branch into the LSA vessel via the ostium. Push the handle to advance the Delivery System until the collar of the Arcevo™ LSA implant aligns with the transected thoracic aorta, the LSA branch is fully inserted within the LSA vessel, and the implant is held against the outer curve as shown in **Figure 10**.

CAUTION: DO NOT manipulate the Arcevo™ LSA system excessively when aligning the collar to the transected aorta. Excessive manipulation or advancing the device further into the aorta may result in LSA branch migration and loss of seal of the LSA branch within the LSA vessel. Ensure the LSA branch is positioned within the LSA vessel.

11.3.4 Fully insert the LSA branch into the LSA vessel by slowly pushing the handle in the direction of the outer curve. If excessive resistance is felt when introducing the Arcevo™ LSA delivery system with loaded implant, DO NOT advance further. Stop and assess the source of the resistance before proceeding further.

WARNING: DO NOT continue to advance the device during the insertion of the Arcevo™ LSA delivery system if excessive resistance is felt, as the risk of vessel trauma or delivery system damage may occur. Additional care must be taken in case of acute dissection to avoid entering the false lumen. Stop and assess the cause of resistance before proceeding with the Arcevo™ LSA procedure.

11.4 STENT GRAFT DEPLOYMENT AND DELIVERY SYSTEM WITHDRAWAL

11.4.1 Ensure that the Arcevo™ LSA implant is fully inserted into the anatomy before initiating deployment (e.g., the Arcevo™ LSA collar should be aligned with the transected aorta and the LSA branch should be fully inserted into the LSA vessel).

11.4.2 Hold the handle steady and ensure the stented LSA branch and suture collar maintain position throughout deployment. While deploying the aortic and LSA release mechanisms it is recommended to secure the Arcevo™ LSA implant by gently holding the device and the aortic wall in place at the level of the transected aorta with forceps. This ensures the device stays in the intended location during deployment and delivery system withdrawal.

11.4.3 Deploy the aortic stent graft by first pulling the trigger (**Figure 11A**) and then retracting the handle (**Figure 11B**) until the delivery system locks to open the Wrap and deploy the stent graft (**Figure 12**).


Figure 11A. Pull aortic trigger


Figure 11B. Retract handle until hard stop

11.4.4 Deploy the LSA branch second by pulling the LSA Pull Tab (**Figure 13**) which will unlace the LSA Wrap and deploy the stented branch.

CAUTION: DO NOT completely remove the LSA Deployment Wire from the delivery system; removing the wire may dislodge the implant or result in delivery system damage.

11.4.5 Continue to retract the LSA Deployment Wire via LSA Pull Tab until resistance is felt and the LSA Tip retracts toward the Handle to assist with Delivery System withdrawal.


Figure 12. Deployment of Arcevo™ LSA implant in transected aorta after positioning


Figure 13. Retract the LSA Pull Tab to deploy the LSA branch after aortic deployment

11.4.6 Once the Arcevo™ LSA implant is fully deployed, remove the Shaping Wire, if an aortic guidewire was not used.

11.4.7 The Arcevo™ LSA delivery system can now be carefully withdrawn with slight rotational movements, if needed.

11.5 PROCEDURE COMPLETION –ANASTOMOSES

11.5.1 Once the Arcevo™ LSA delivery system has been removed, the collar must be sutured to the native aorta. Complete the anastomosis of the Arcevo™ LSA collar to the transected aorta using sterile surgical suture and standard surgical techniques (**Figure 14**).

CAUTION: Failure to seal the Arcevo™ LSA implant by complete anastomosis of the collar may result in bleeding and/or treatment failure.

11.5.2 Implant a surgical graft in the ascending aorta per institutional protocol and the graft manufacturer's Instructions for Use. Perform the anastomosis of a surgical graft (not supplied) to the transected aorta and Arcevo™ LSA implant using standard surgical techniques. Discretionary branch vessel anastomoses can be completed in standard surgical fashion per institutional protocol. See **Figure 15** for representative depiction of completed procedure of Arcevo™ LSA with distal landing zone.


Figure 14. Distal anastomosis of Arcevo™ LSA implant by suturing collar to thoracic aorta


Figure 15. Representative total arch replacement with Arcevo™ LSA implant

12.0 DISTAL PLACEMENT OF A THORACIC EXTENSION STENT GRAFT

If the Arcevo™ LSA implant requires extension to complete treatment of the aneurysm or dissection, the Medtronic Valiant™ thoracic stent graft (Closed Web straight configuration or Closed Web tapered configuration only) with the Captivia™ delivery system should be used in a future endovascular procedure based on available in-vitro compatibility testing of both devices under clinically relevant conditions and femoral approach. Implantation of the Medtronic Valiant™ thoracic stent graft shall be performed in accordance with the manufacturer's Instructions for Use. During the distal extension endovascular procedure, take care to manage the position of the Arcevo™ LSA implant, as the insertion and positioning of the Medtronic Valiant™ thoracic stent graft may potentially displace or dislodge the Arcevo™ LSA implant. Observation of the position, alignment, and overlap region of both devices is required under appropriate imaging (e.g., fluoroscopy).

CAUTION: When using the Medtronic Valiant™ thoracic stent graft, the minimum recommended length of overlap between devices is 50 mm as indicated by the chain shaped radiopaque marker. Less than 50 mm of overlap length may result in endoleak or component separation. Avoid excessive overlap in the aortic arch.

The minimum recommended length of overlap between the Arcevo™ LSA implant and the Medtronic Valiant™ thoracic stent graft is 50 mm as indicated by the chain shaped radiopaque marker. For overlap zones in a dissection, adhere to the aortic sizing guidelines in the Medtronic Valiant™ Thoracic Stent Graft with the Captivia™ Delivery System Instructions for Use (but not larger than 10% oversizing compared to the Arcevo™ LSA stent graft nominal aortic diameter).

In patients with aneurysm, the distal extension procedure should be planned and conducted to ensure that both Arcevo™ LSA and Medtronic Valiant™ Captivia™ thoracic stent grafts take the outer curve of the aorta between the proximal anastomosis (above the aneurysm) and a healthy section of the distal thoracic aorta (below the aneurysm). When the distal end of the Arcevo™ LSA is implanted unsupported within the aneurysm sac (with no distal seal until the thoracic stent graft is placed), there may be an increased risk of thrombus generation until the Medtronic Valiant™ thoracic stent graft is implanted. For these patients it is recommended that a Medtronic Valiant™ thoracic stent graft with a proximal diameter that is 4-5 mm greater than the nominal distal diameter of the Arcevo™ LSA be used.


Figure 16. Image of Arcevo™ LSA ring radiopaque markers (solid arrows) and overlap radiopaque marker (open arrow) for distal extension using Medtronic Valiant™ Captivia™ thoracic stent graft

13.0 PATIENT FOLLOW-UP

Note: The follow-up recommendations below apply to the Arcevo™ LSA Hybrid Stent Graft System only. If the Medtronic Valiant™ Captivia™ thoracic stent graft is used to extend the Arcevo™ LSA implant distally, refer to the Medtronic Valiant™ Thoracic Stent Graft with the Captivia™ Delivery System Instructions for Use for the recommended follow-up procedure.

All patients should be counseled that a total arch repair procedure with the Arcevo™ LSA Hybrid Stent Graft System requires lifelong periodic follow-up to assess their health, along with imaging to assess implant performance, as advised by their physician. Patients should be advised on the importance of adhering to the follow-up schedule within the first year of their procedure and annually thereafter, to ensure ongoing safety and effectiveness of their aneurysm and/or dissection repair. Ultimately, it is the physician's responsibility to determine the appropriate imaging schedule for an individual patient based on previous clinical outcomes and the overall patient clinical condition.

Current recommended imaging modalities for patients with an implanted stent graft include X-Ray and Computed Tomography (CT) with or without contrast medium; alternate modalities such as magnetic resonance imaging (MRI)/angiography (MRA) duplex ultrasound may be considered in patients with impaired renal function or those allergic to contrast medium, at physician discretion.

Post-operative imaging is recommended to include a non-contrast CT, followed by an early and late phase contrast CT (slice thickness < 3 mm recommended) to provide information on endoleaks, stent graft migration, component separation, thrombosis, and morphological changes. CT imaging is recommended to be performed at predischARGE or at 1 month, 3 - 6 months, 12 months and yearly thereafter, if stable condition. Magnetic Resonance Angiography (MRA) may be used to avoid cumulative lifetime radiation exposure and iodinated contrast media. After three years follow-up, the interval of CT or MRI imaging can be extended to two (2) to three (3) years in patients with stable aneurysm who have no endoleak and no other sign of complications.

Patients observed with specific clinical conditions, such as enlarging aneurysms or false lumen (e.g., aortic dilation > 5 mm with or without endoleak), changes in the structure, seal, or position of their stent graft, should receive additional follow-up. Patients with early or late type I and III endoleak, enlarging aneurysm, or stent graft migration should undergo prompt re-intervention. Endovascular procedures should be taken into consideration. Patients with type II or type IV endoleak or aortic growth > 10 mm should receive enhanced follow-up with more frequent intervals. Additional CTA (MRA) are recommended at three (3) and six (6) months after detection. Care should be taken during planning of follow-up in younger patients regarding the cumulative lifetime radiation exposure caused by multiple CT examinations. A combination of multi-view chest X-rays and MRA may be an alternative long-term follow-up in young patients.

14.0 SUMMARY OF CLINICAL STUDIES AND CLINICAL BENEFIT

The Arcevo™ LSA Hybrid Stent Graft System is an investigational device in the United States and Europe. The outcomes of the clinical study for this device will be provided in a future revision to this IFU.

15.0 REPORTING OF SERIOUS INCIDENTS AND ADVERSE EVENTS

Serious incidents that occur during the Arcevo™ LSA procedure or after Arcevo™ LSA implantation should be reported per the ARTIZEN study protocol.

16.0 DISCLAIMER OF WARRANTIES; LIMITS OF LIABILITY

EXCEPT AS EXPRESSLY SET FORTH IN WRITING BY ARTIVION, ARTIVION MAKES NO REPRESENTATIONS OR WARRANTIES OF ANY KIND, EXPRESS OR IMPLIED, WITH RESPECT TO THIS PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL ARTIVION BE LIABLE FOR ANY MEDICAL EXPENSES, OR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, MISUSE, ALTERATION OR UNAUTHORIZED REPAIR, DEFECT, FAILURE, OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT, OR OTHERWISE. In the event that any part or term of this disclaimer is found invalid or unenforceable for any reason: (i) the validity of the remaining portions of this disclaimer and limitation shall not be affected; (ii) any action for breach of warranty must be commenced within one year after any such claim or cause of action accrued; and (iii) the remedy for any such breach is limited to the replacement of the product.

17.0 SYMBOLS USED ON LABELING

	Manufacturer		Use By Date		Keep Dry		Patient Identification
	Date of Manufacture		Consult Instructions for Use		Temperature Limit		Date of Implantation
	Catalog Number		Sterilized Using Ethylene Oxide		MR Conditional		Name and Address of the implanting healthcare institution/ provider
	Model Number		Do Not Reuse		European Authorized Representative		Caution
	Lot Number		Do Not Re-sterilize		Information Website for Patients		Single Sterile Barrier System
	UDI as AIDC format		Do Not Use If Package Is Damaged		Medical Device		Double Sterile Barrier System

18.0 PATIENT CARD

The Arcevo™ LSA Hybrid Stent Graft System is packaged with a Patient Implant Card, all patients should be instructed to keep this card in their possession at all times.

Instructions for completion of Patient Implant Card (to be completed by the healthcare institution/provider):

1. Patient Identification: Patient name (first, middle, last) and/or patient ID
2. Date of implantation (YYYY-MM-DD)
3. Name and address of the implanting healthcare institution/provider and/or Physician's name

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
Telefon: +1 888.427.9654 • Fax: 770.590.3753 • Email: customerservice.us@artivion.com

Artivion, Inc.
1655 Roberts Blvd. NW
Kennesaw, GA 30144 USA
+1 (888) 427-9654
www.artivion.com

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