


## Ordering Information

Diameter of Surgical Wound Opening	CO <sub>2</sub> Flow Recommended	Recommended Product	Catalog Name	Contents	Reference Image
<3 cm	3 L/Min	CarbonMini CO <sub>2</sub> Diffuser	CARBONMINI	8 Units / Box	
3-4 cm	5 L/Min				
5-8 cm	8 L/Min				
>8 cm	10 L/Min	CarbonAid CO <sub>2</sub> Diffuser	CARBONAID	8 Units / Box	