

Ordering Information

Chord-X Pre-Measured Loops and Adjustable Suture System

ePTFE Suture Diameter: USP 2-0 (Diameter similar to GORE-TEX CV-4)	18mm Taper Point Needle Options, Catalog Numbers	
Loop Length	3/8 Circle	1/2 Circle
Adjustable	CXL-20-1838-0	CXL-20-1812-0
12mm	CXL-20-1838-12	CXL-20-1812-12
16mm	CXL-20-1838-16	CXL-20-1812-16
20mm	CXL-20-1838-20	CXL-20-1812-20
24mm	CXL-20-1838-24	CXL-20-1812-24

- 1 box = 5 sterile Chord-X packaged units.
- Each unit contains a Chord-X prosthesis featuring one suture pair for the papillary muscle and three independent pre-measured chordal suture loops or suture strands for leaflet attachment.
- Each prosthesis is positioned on disposable foam base-plates for suture management.
- Each unit has double-armed taper point needles, two PTFE pledgets: one attached and one independent.

Chord-X Chordal Sizer

Packaged Quantity	Use	Catalog Number
1 box = 5 individually packaged sterile Chord-X Chordal Sizers	Single use, disposable device	CXCS