

On-X[®] Ascending Aortic Prosthesis Configurations

Valve Size	Catalog Number	Valve Internal Diameter (mm)	Sewing Ring Diameter (mm)	Graft Internal Diameter (mm)	Graft Length (cm)
19	ONXAAP-19	17.4	27	20	11
21	ONXAAP-21	19.4	30	22	11
23	ONXAAP-23	21.4	33	24	11
25	ONXAAP-25	23.4	34	26	11
27/29	ONXAAP-27/29	23.4	36	26	11

ARTIVION™